

# Dental Instruments Catalog

THEY USE



Manufacturers & Exporters of Quality Dental & Surgical instruments

**WE MAKE TOOLS**



### Dear Valued Buyers

We feel immense pleasure to present our latest catalogue in which we have featured a most comprehensive range of all types of quality Surgical & Dental instruments required for day to day use for Hospital. We will be pleased to quote you a competitive Price and Delivery time.

All instruments are manufactured according to the recognized standard patterns currently in use in Hospital throughout the world and are being produced under restrict control of experts to avoid any flaw and to refrain any complaint at any stage.

A group of very vast experienced persons in manufacturing handling our quality Control System.

Any instruments proven defective in workmanship or material will be replaced with out charges.

We express our thanks to all our customers for the confidence they have placed us in the past. We further assure them that we will continue our best efforts to gain their confidence in the future.

We are always at your service.

Chief Executive



# CERTIFICATE OF REGISTRATION



## Certificate of Registration



The Governing Board of  
Q.A. International Certification Limited  
hereby grants to:

**AZQA INSTRUMENTS**  
Registration No: QAIK / PK / 1712 - B

hereafter called the Registered Company the right to be listed in the Directory of Registered Companies in respect of the services listed below. These services shall be offered by the Registered Company at or from only the address given below in accordance with the quality management system in compliance with the Requirements of ISO 13485:2016.

Address to which this Certificate refers:

Sadey Shehar Road, SIE Sialkot 51310 Pakistan

Approved Scope to which this Certificate refers:

**Manufacturers of Non-Acute Surgical, Dental Instruments & Scissors.**  
(Please note that the above scope represents the certified activity of the named organization and so each organization may undertake additional activities that are not covered under this certification).

Signed for and on behalf of the Board

*[Signature]*  
CHIEF EXECUTIVE

*[Signature]*  
SCHEME MANAGER

Certificate Issue Date: 15th February 2019 - Certificate Renewal Before: 15th December 2019  
Date of Initial Registration: 15th December 2017 - Re-Certification Before: 15th December 2020

This Certificate of Registration is granted subject to the Regulations approved by the Board.

QA INTERNATIONAL

UK Accredited Certification Ltd  
Buckley Court  
Buckley Road  
Buckley  
Wrexham, Cheshire  
WA9 7JL, UK  
Tel: +44 (0)1246 594073  
Fax: +44 (0)1246 594080  
www.qai.org.uk



## Certificate of Registration



The Governing Board of  
Q.A. International Certification Limited  
hereby grants to:

**AZQA INSTRUMENTS**  
Registration No: QAIK / PK / 1813 - A

hereafter called the Registered Company the right to be listed in the Directory of Registered Companies in respect of the services listed below. These services shall be offered by the Registered Company at or from only the address given below in accordance with the quality management system in compliance with ISO 9001:2015.

Address to which this Certificate refers:

Sadey Shehar Road, SIE Sialkot 51310 Pakistan

Approved Scope to which this Certificate refers:

**Manufacturers of Surgical, Dental, Veterinary, Beauty Care Instrument & Scissors.**  
(Please note that the above scope represents the certified activity of the named organization and so each organization may undertake additional activities that are not covered under this certification).

Signed for and on behalf of the Board

*[Signature]*  
CHIEF EXECUTIVE

*[Signature]*  
SCHEME MANAGER

Certificate Issue Date: 25th February 2019 - Certificate Renewal Before: 15th February 2020  
Date of Initial Registration: 15th February 2019 - Re-Certification Before: 15th February 2022

This Certificate of Registration is granted subject to the Regulations approved by the Board.

QA INTERNATIONAL

UK Accredited Certification Ltd  
Buckley Court  
Buckley Road  
Buckley  
Wrexham, Cheshire  
WA9 7JL, UK  
Tel: +44 (0)1246 594073  
Fax: +44 (0)1246 594080  
www.qai.org.uk



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

APR 16 2019

Mr. Muhammad Riaz  
Chief Executive  
AZQA INSTRUMENTS  
Sadey Shehar Road,  
SIE Sialkot  
51310 Pakistan

and

Mr. Muhammad Ashraf  
Quality Management Consultant  
Near Model Masjid  
Nasir Road  
Sialkot, Pakistan

Dear Messrs. Riaz and Ashraf:

This is to acknowledge receipt of a March 1, 2019, letter from Mr. Muhammad Ashraf certifying the compliance of AZQA INSTRUMENTS Industries (AZQA) with the Food and Drug Administration (FDA) Quality System Regulation of 1997, which includes the current good manufacturing practice (CGMP) requirements. The Quality System Regulation is set forth in Title 21, Code of Federal Regulations (CFR), Part 820. The consultant certification confirmed that a quality system audit of Maymoy was performed February 28, 2019.

The quality system audit report states that manufactures reasonable surgical instruments. Based on our review of the audit results and certification, (AZQA) has been placed on Attachment A of Import Alert 876-01. You may begin exporting devices to the United States (U.S.) that were manufactured after the consultant certified your firm's compliance with the CGMP's; however, your shipments may be subject to the guidance outlined in Attachment A of Import Alert 876-01. After five consecutive shipments comply with the import alert guidance, you may request your firm be placed on Attachment B. Submit your request directly to the FDA district office for their concurrence and further submission to this office for action.

The placement of the firm on Attachment A is limited to devices manufactured under the name of AZQA INSTRUMENTS Sadey Shehar Road, SIE Sialkot 51310 Pakistan.

In the event the manufacturing name and/or address change, FDA requests that notification be immediately forwarded to this office. A change in the name and/or address of the manufacturing facility without notifying FDA will result in a re-evaluation of the compliance status of your firm.

Page 2 - Messrs. Riaz and Ashraf

The FDA is assured by both the consultant and (AZQA) that no devices will be imported into the U.S. that are not specifically manufactured by AZQA. This includes other manufacturers in Pakistan completing specific manufacturing processes for AZQA, and other manufacturers in Pakistan providing their devices to AZQA for importation into the U.S.

The decision based on your consultant certification will remain in effect until such time as FDA is able to visit Sialkot, Pakistan for an inspection of your facility. During this inspection all corrections and procedures will be evaluated and confirmed. Any new CGMP deviations, or any uncorrected deviations that were previously certified to, may result in a re-evaluation of the compliance status of your firm, AZQA, including the possibility of removal from Attachment A.

We request that a quality system follow up audit be performed at AZQA, within six months of exporting devices to the U.S. You will be advised of the timing of FDA's inspection schedule.

AZQA has an ongoing responsibility to conduct internal self-audits to assure you continue to maintain conformance with the Quality System Regulation.

If you have any questions regarding this correspondence, or need further assistance, please contact me at (301) 594-4595 or FAX (301) 594-4636.

Sincerely yours,

*[Signature]*  
Carol J Shirk

Carol J Shirk  
Consumer Safety Officer  
General Surgery Devices Branch  
Division of Enforcement I  
Office of Compliance  
Center for Devices and  
Biological Health



Medicines and Healthcare products Regulatory Agency

Medicines and Healthcare products  
Regulatory Agency  
Milton Keynes  
1 Nine Elms Lane London W8 5NF (UK)  
General enquiries  
Telephone: 020 7044 5100 Fax: 020 7044 5141  
Email: [ask@mhra.gov.uk](mailto:ask@mhra.gov.uk)  
www.mhra.gov.uk  
Helpline: 0800 135 0810  
Direct Fax: 0207 044 1442  
Email: [web@mhra.gov.uk](mailto:web@mhra.gov.uk)  
Charity: CA 200803

Mr D Iqbal  
Blue Star Company  
48, Therncliffe Road  
Newstead Green  
Southall  
Middlesex  
UB8 3RQ

09 October 2016

Dear Mr. D Iqbal,

MEDICAL DEVICES REGULATION 2002: REGULATION 19  
Registration of Persons Placing Medical Devices on the Market

Thank you for informing the Competent Authority of the details of Manufacturers Name:-  
located at Manufacturers Address:-  
for whom you are acting as the authorized representative and for supplying the medical device information.

Your registration has been recorded based on your declaration that you have determined that the device(s) fall within the definition of "medical device", and that you have classified it/them as falling within Regulation 19 taking into account the intended purpose(s) and mode(s) of action. In accepting your registration, I should make clear that the Competent Authority does not examine each individual notification and therefore cannot and does not necessarily endorse these determinations. Neither does this letter represent any form of accreditation or approval by the UK Competent Authority.

Your registration is based upon your declaration on the RG2 form and means that:

**For Manufacturers of Class I medical devices, Assemblers, and Distributors**  
You should now be operating under the Medical Device Directive and the above Regulations for the products you asked us to register by fully complying with the essential requirements. CE marking these products or labelling them as such.

**For Manufacturers of Custom-made devices**  
You should be ready to obtain compliance with the Directive and Regulations and should be manufacturing custom-made devices in accordance with their requirements.

If you stop placing devices on the market or if you are not complying with the Regulations you should inform us so that we can amend our records. You should be aware that it is an offence to place on the market CE marked devices that do not comply with the regulations.



An Executive Agency of the Department of Health

The information you provided has been recorded against the reference number shown at the top of this letter, which we ask you to quote in all future correspondence and communications.

Please inform us of any changes to:

- The company information
- Additional generic group of devices (any individual products within an existing generic group)
- Discrimination of a generic group of devices

Please use RG2, the Registration form, to tell us about any of these changes.

Thank you for registering the following generic groups of devices:

**Class I Devices:**  
Dental Instruments (Re-Usable & Non-Powered)  
Surgical Instruments (Re-Usable And Non-Powered)

**Custom Made Devices:**

None

**Products Covered By Article 12:**

None

Should you have any queries regarding your registration please contact us on the telephone number given at the top of this letter.

Yours sincerely

*[Signature]*  
Rob Higgins



CA014167



ACCREDITED COMPANY  
BS EN ISO 9001:2008  
BS EN ISO 13485:2003  
QIA/CPK/3587-A  
QIA/CPK/3588-B



300957719



GMP Approved Facility



A-20939



S286



1256696-9

All third party names/logos/trademarks are the property of their respected owners.

E-mail: [azqa.instruments@gmail.com](mailto:azqa.instruments@gmail.com)



[info@azgainstruments.com](mailto:info@azgainstruments.com)

**01****EXTRACTING FORCEPS****Page # 01 - 28****02****ROOT ELEVATORS****Page # 29 - 48****03****DIAGNOSTIC, PROBES, SCALERS, FILLING  
& LABORATORY INSTRUMENTS****Page # 49 - 98****04****MAXILLOFACIAL SURGERY  
INSTRUMENTS****Page # 99 - 124****05****ORTHODONTIC INSTRUMENTS  
& PLIERS****Page # 125 - 144****06****DENTAL IMPLANT INSTRUMENTS  
& Impression Tray****Page # 145 - 182****07****HOLLOWARE****Page # 183 - 191**



## Handle Selection For Extracting Forceps

Dear Valued Buyers if any Universal Fig. No. of Extracting Forceps not Found in our Catalogue you Can Refer the Fig. No. From These Diagram.



Dear Valued Buyers these are our Handle Styles



English Pattern



English Pattern



English Pattern, Anatomical Handle



American Pattern



Deep Grip Forceps With Hole Pattern



Hole Pattern



## EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)

Function: To remove teeth from bony socket



**DE-001**  
Fig. 1  
Upper Centrals  
and Canines



**DE-002**  
Fig. 2  
Upper Laterals  
and Canines



**DE-003**  
Fig. 4  
Lower Incisors  
and Canines



**DE-004**  
Fig. 7  
Upper Premolars



**DE-005**  
Fig. 8  
Lower Premolars



**DE-006**  
Fig. 13  
Lower Premolars



**DE-007**  
Fig. 17  
Upper Molars,  
Right



**DE-008**  
Fig. 18  
Upper Molars, Left



**DE-009**  
Fig. 21  
Lower Molar



**DE-010**  
Fig. 22  
Lower Molar



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



# EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)

Function: To remove teeth from bony socket



**DE-011**  
**# 18A**  
upper molars



**DE-012**  
**# 19**  
upper third molars



**DE-013**  
**# 20**  
lower third molar



**DE-014**  
**# 24**  
lower molar, right



**DE-015**  
**# 29**  
upper roots



**DE-016**  
**# 30**  
Upper Roots



**DE-017**  
**# 33**  
lower roots



**DE-018**  
**# 33A**  
Lower Roots



**DE-019**  
**# 33M**  
Lower roots  
with long beak



**DE-020**  
**# 46L**



**DE-021**  
**# 51**  
Upper Roots



**DE-022**  
**# 51A**  
Upper Roots



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines

## EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)

Function: To remove teeth from bony socket

### WISDOM EXTRACTING FORCEPS



**DE-023**  
**Routerier**  
**(# 22 1/2 L)**  
lower molars and  
wisdoms, left



**DE-024**  
**Routerier**  
**(# 22 1/2 R)**  
lower molars and  
wisdoms, right



**DE-025**  
**Routerier**  
**(# 33 1/2 L)**  
lower molar and  
wisdom roots, left



**DE-026**  
**Routerier**  
**(# 33 1/2 R)**  
lower molar and  
wisdom roots, right



*Routerier extracting forceps for lower molars  
and third molars, right*



**DE-027**  
**Fig. 186L**  
lower molars



**DE-028**  
**Fig. 186R**  
lower molars



**DE-029**  
**# 51AL**



**DE-030**  
**# 67**  
Upper Third Molars



**DE-031**  
**# 67A**  
Upper Third Molars



**DE-032**  
**# 79**  
Lower Third Molars



*Angled extracting forceps for upper roots*




## EXTRACTING FORCEPS (English Pattern)


- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)



**DE-033**  
# 73  
Lower Molars



**DE-034**  
# 73 A  
lower molars



**DE-035**  
# 73 S



**DE-036**  
# 74  
Lower Roots



**DE-037**  
# 74N  
Lower Roots



**DE-038**  
# 75  
Lower Premolars



**DE-039**  
# 86A  
Lower Molars



**DE-040**  
# 86B  
lower molars



**DE-041**  
# 86C  
lower molars



**DE-042**  
# 87  
Lower Molars



**DE-043**  
# 89  
Upper Molars,  
Right



**DE-044**  
# 90  
Upper Molars,  
Left



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines

## EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)

### A-TRAUMATIC EXTRACTION

Improved tapered beaks for complete root contact. Cross serrations. Deep Gripe Extracting Forceps



• No Fracture

• No Tissue Damage

• Extraction With Easy

Thin Tip



**DE-045  
# 34**

Upper incisors and  
premolars  
deep-gripping

Wide Tip



**DE-046  
# 34A**

Upper centrals / cains  
deep-gripping

Thin Tip



**DE-047  
Fig. 35**

upper permolars  
deep-gripping

Wide Tip



**DE-048  
# 35A**

Upper incisors and  
premolars  
deep-gripping

Thin Tip



**DE-049  
# 36**

Lower premolars  
deep-gripping

Wide Tip



**DE-050  
# 36A**

Lower premolars  
deep-gripping



**DE-051  
# 336**

lower permolars  
deep-gripping



**DE-052  
F-7C**

lower molar



**DE-053  
# 9**

Apex Forceps Lower Universal  
Deep Root Pick



**DE-054  
# 79N**

Lower third molars  
deep-gripping



**DE-055  
# 79**

Lower Thrid Molars



#### Function:

Extraction forceps are used to further loosen teeth by rotational and tilting movements and finally to extract the teeth in combination with a fulcrum.



## EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Давиеры (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)

### DEEP ROOT PICK FORCEPS



**DE-056**  
**# 146**

Apex Forceps Lower Root  
Deep Root Pick



**DE-057**  
**# 153**

Apex Forceps Upper Root  
Deep Root Pick



**DE-058**  
**# 149**

Apex Forceps Upper Root  
Deep Root Pick



**DE-061**  
**Fig. 44**



**DE-060**  
**Fig. 44Z**



**DE-059**  
**Fig. 45Z**



**DE-062**  
**# 46XL**  
Apex Forceps Lower Root  
Deep Root Pick



**DE-063**  
**# 67**  
Upper Third Molars



**DE-064**  
**# 67A**  
Upper Third Molars



forceps for upper incisors, premolars and roots



**DE-065**  
**# 67N**  
Lower third molars  
deep-gripping



**DE-066**  
**# 331**  
upper roots  
deep-gripping

## EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



Fig. 1  
upper  
laterals  
and canines



**DE-067**  
Set of 10 Pieces  
Fig. 1  
Fig. 7  
Fig. 13  
Fig. 17  
Fig. 18  
Fig. 22  
Fig. 51A  
Fig. 67  
Fig. 74N  
Fig. 79

Fig. 7  
upper premolars



Fig. 13  
lower premolars



Fig. 17  
upper molars, right



Fig. 18  
upper molars, left



Fig. 22  
lower roots



Fig. 51A  
upper roots



Fig. 67  
upper third molars



Fig. 74N  
lower molar



Fig. 79  
lower third molars



Extracting forceps for upper roots





- Щипцы экстракционные
- Daviers
- Zahnzangen

- Щипці екстракційні
- Forceps de Extracción
- استخراج الملقط

## Pediatric Tooth Extracting Forceps Set



**DE-069**  
English Pattern Klein #51S  
Pedodontic Upper root



**DE-068**  
English Pattern Klein #137  
Pedodontic, Upper Incisors Spring handle



**DE-072**  
English Pattern Klein #5  
Pedodontic Lower Incisors.



**DE-070**  
English Pattern Klein #139  
Pedodontic Upper Premolars

### Baby Tooth Extraction Forceps




**DE-075**  
Fig. # 01 to 07  
Set of 7 Pieces



**DE-073**  
English Pattern Klein #7  
Pedodontic Lower Roots,



**DE-071**  
English Pattern Klein #3  
Pedodontic Upper Molars Universal



**DE-074**  
English Pattern Klein #6  
Pedodontic Lower Molars, Spring Handle



Witzel  
**DE-076**  
170mm  
Upper roots, Universal



Witzel  
**DE-077**  
165mm  
Lower roots, Universal



**DE-078**  
#400(300)  
165mm  
Upper roots, Universal



**DE-079**  
#401(301)  
160mm  
Lower roots, Universal



**DE-080**  
Witzel Forceps With Gold  
Set of 3 Pcs



**DE-084**  
Witzel Forceps With Hole And Screw Pattern  
Set of 3 Pcs



**DE-088**  
Witzel Forceps With Anatomical Handle  
Set of 3 Pcs



**DE-081**  
#401(301)  
160mm  
Lower roots, Universal



**DE-085**  
#401(301)  
160mm  
Lower roots, Universal



**DE-089**  
#401(301)  
160mm  
Lower roots, Universal



**DE-082**  
#400(300)  
165mm  
Upper roots, Universal



**DE-086**  
#400(300)  
165mm  
Upper roots, Universal



**DE-090**  
#400(300)  
165mm  
Upper roots, Universal



Witzel  
**DE-083**  
165mm  
Lower roots, Universal



Witzel  
**DE-087**  
165mm  
Lower roots, Universal



Witzel  
**DE-091**  
165mm  
Lower roots, Universal



## ROOTS SPLINTER FORCEPS

- Щипцы для удаления корневых осколков
- Pincers Fragments
- Wurzelsplitterzangen

- Щипці для видалення корневих осколків
- Pinzas Para Fragmentos
- جنور منشقة الملقط



**DE-092**  
Splinter Forceps  
15.5cm



*Stieglitz Root Splinter Forceps*




**DE-093**  
Splinter Forceps  
15.5cm 45°



**DE-094**  
Stieglitz  
Root Fragments




**DE-095**  
Stieglitz  
Root Fragments



**DE-096**  
Stieglitz Splinter Forceps  
140mm 90°



**DE-097**  
Stieglitz Splinter Forceps  
140mm Straight




**DE-098**  
130mm

With Diamond



**DE-099**  
130mm

With Diamond




**DE-100**  
125mm


With Diamond



**DE-101**  
130mm



**DE-102**  
130mm



**DE-103**  
125mm



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Dавiers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- ستخراج الملقط (النمط الأمريكي)



**DE-104**  
**# 1STD**  
Winter  
upper centrals and canines

**DE-105**  
Fig. 1  
upper incisors  
and roots

**DE-106**  
Cryer  
Fig. 150

## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzeper per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المنقط (النمط الأمريكي)





**DE-108**  
Fig. 62  
incisors, premolars,  
deciduous teeth



**DE-109**  
# 10S  
upper roots



**DE-110**  
# 17  
lower molars



**DE-111**  
Harris  
Fig. 18R



**DE-112**  
Harris  
Fig. 18L




**DE-107**  
Fig. 150A  
Cryer  
upper incisors, premolars, roots



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- ستخراج الملقط (النمط الأمريكي)



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- استخراج المنقط (النمط الأمريكي)



**DE-117**  
Fig. 24  
upper molars



**DE-118**  
# 210S  
upper third molars



**DE-119**  
# 210H  
upper third molars



**DE-120**  
# 286  
Upper Roots,  
Premolars, Incisors



**DE-123**  
# 286  
Lower Molar



**DE-121**  
# 99C  
Kells  
upper incisors,  
canines, premolars



**DE-122**  
# 222  
lower third molars



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Dавiers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- ستخراج الملقط (النمط الأمريكي)



**DE-124**  
Nevius  
Fig. 88L  
upper molars, left



**DE-125**  
Fig. 65  
upper incisors  
and roots



**DE-126**  
Parmly Fig. 32A  
upper premolars  
and roots




**DE-127**  
Parmly Fig. 32  
upper premolars  
and molars



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Daviers - forme américaine
- Pinzeper per estrazione - forma americana
- Fórceps de extracción - forma americana
- استخراج المنقط (النمط الأمريكي)



**DE-128**  
COHEN  
universal

**DE-129**  
Fig. 203  
lower incisors  
and roots

**DE-130**  
Tomes Fig. 69  
upper and lower  
roots



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Давиеры - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- استخراج الملقط (النمط الأمريكي)



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- استخراج المنقط (النمط الأمريكي)



**DE-136**  
Hull  
Fig. 101  
Deciduous teeth

**DE-137**  
Cryer  
Fig. 150S  
upper teeth,  
for children

**DE-138**  
Fig. 150 AS

## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- ستخراج الملقط (النمط الأمريكي)



**DE-139**  
Cryer  
Fig. 151  
lower teeth,  
for children

**DE-140**  
Fig. 151 S  
lower incisors,  
premolars, roots

**DE-141**  
Fig. 151 AS



## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzper per estrazione - forma americana
- استخراج المنقط (النمط الأمريكي)




**DE-142  
# 16S**  
lower molars for children



**DE-143  
# 17S**  
lower molars



**DE-144  
# 23 S**  
Lower molars



**DE-145  
# 27**  
lower molars



**DE-146  
# 18LS**  
upper molars left



**DE-147  
# 18RS**  
upper molars right



**DE-148  
# 85 Tomas**  
lower canines, premolars and molars



**DE-149  
# 85A Tomas**  
lower canines, premolars and molars



**DE-150  
# 69S**  
upper and lower roots



**DE-151  
# 99A Kells**  
upper incisors, canines, premolars



**DE-152  
# 217**  
lower molars

## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Dавiers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- ستخراج الملقط (النمط الأمريكي)



**DE-153**  
**#10 SK**  
Upper Molars Sup



**DE-154**  
**#17 SK Universal**  
Lower Molars



**DE-155**  
**#23 SK Universal**  
Upper Molars



**DE-156**  
**#69 SK Tomas**  
upper and lower roots



**DE-157**  
**#99 SK**  
upper incisors,  
canines, premolars



**DE-158**  
**#150 SK**  
upper incisors,  
premolars, roots



**DE-159**  
**#151 SK**  
lower incisors,  
premolars, roots



**DE-160**  
**#222 SK**  
lower third molars

## EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- استخراج المنقط (النمط الأمريكي)



**DE-161**  
Fig. 1  
Upper Anterior

**DE-162**  
Lower Molar  
Fig. 23 Cow Horn



**DE-163**  
Fig. 74  
Lower Anterior



**DE-164**  
Fig. 88L  
Upper Molar Left



**DE-165**  
Fig. 88R  
Upper Molar Right



**DE-166**  
Fig. 74  
Lower Anterior



**DE-167**  
Fig. 106  
Lower Incisors



**DE-168**  
Fig. 174  
Extra Narrow  
Lower Anterior  
Apical Retention  
Forceps



**DE-169**  
Fig. 859  
Lower Molars  
Roots





# EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Щипці екстракційні (американська форма)
- Pinzeper per estrazione - forma americana
- Férceps de extracción - forma americana
- استخراج الملقط (النمط الأمريكي)

## Carbide Diamond Coated Point



**DE-170**  
150 T.C.  
Carbide Serrated Beak  
Upper Incisors, Cuspids,  
Bicusps & Roots,  
Universal



**DE-171**  
151 T.C.  
Carbide Serrated Beak  
Lower Incisors, Cuspids,  
Bicusps & Roots, Universal



**DE-172**  
1STD T.C.  
Carbide Serrated Beak  
Upper Incisors, Cuspids,



**DE-173**  
17 T.C.  
Carbide Serrated Beak  
1st & 2nd Lower Molars,  
Universal



**DE-174**  
53L T.C.  
Carbide Serrated Beak  
1st & 2nd Upper Molars,  
Left



**DE-175**  
53R T.C.  
Carbide Serrated Beak  
1st & 2nd Upper Molars,  
Right



**DE-176**  
65 T.C.  
Carbide Serrated Beak  
Upper Incisors, Roots &  
Fragments



**DE-177**  
69 T.C.  
Carbide Serrated Beak  
Upper & Lower Fragments  
And Roots



**DE-178**  
150S T.C. Pedro forcep  
Carbide Serrated Beak  
Pedo, Upper Primary Teeth  
and Roots, Universal



**DE-179**  
151S T.C. Pedro forcep  
Carbide Serrated Beak  
Pedo, Lower Primary Teeth  
and Roots, Universal

## BONE ROUNGER

- Костные кусачки
- Hohlmeißelzangen
- Pincers-gouges
- Кусачки для кісток
- Pinzas-gubias
- روتجورس العظام



**DE-180**  
Blumenthal  
15cm



**DE-181**  
Blumenthal  
15cm



**DE-182**  
Blumenthal  
15cm



**DE-183**  
Friedman  
14cm



**DE-184**  
Friedman  
14.5cm



**DE-185**  
Micro-Friedman  
14cm



**DE-186**  
Micro-Friedman  
15cm



**DE-187**  
Micro-Friedman  
16cm



**DE-188**  
Kleinert-Kutz  
18.5cm



**DE-189**  
Beyer  
18.5cm




**DE-190**  
Mead 1A  
18.5cm



**DE-191**  
Luer  
18.5cm

## EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)



**DE-192**  
Meissner  
Extracting Forceps

**DE-193**  
Meissner  
Extracting Forceps

**DE-194**  
Meissner  
Extracting Forceps

**DE-195**  
Comellas Forceps  
Anatomical Handle



## DENTAL SYRINGES

- Шприцы стоматологические карпульные
- Seringues à usage dentaire
- Dentspritzen

- Шприці стоматологічні карпульні
- Jeringas dentales
- المحاقن الأسنان



DE-196 1.8ml  
DE-197 2.2ml



DE-198 1.8ml  
DE-199 2.2ml



DE-200 1.8ml  
DE-201 2.2ml



DE-202 1.8ml  
DE-203 2.2ml



DE-204 1.8ml  
DE-205 2.2ml



DE-206 1.8ml  
DE-207 2.2ml



DE-208  
# 1001  
Intraligamental  
Anesthesia  
Syringe in  
Black Colour



DE-209 1.8ml  
DE-210 2.2ml  
Angled



DE-211  
Aluminium - Dark 1.8 ml

DE-212 Aluminium - Silver 1.8 ml  
DE-213 Aluminium - Orange 1.8 ml  
DE-214 Aluminium - Green 1.8 ml  
DE-215 Aluminium - Blue 1.8 ml



Type C  
DE-216 1.8ml  
DE-217 2.2ml



Type C  
DE-218 1.8ml  
DE-219 2.2ml



2-Ring  
DE-220 1.8ml  
DE-221 2.2ml




Aspirating Syringe  
DE-222 1.8mm  
DE-223 2.2mm




Cartridge Syringe  
DE-224 1.8mm  
DE-225 2.2mm




Cartridge Syringe  
DE-226 1.8mm  
DE-227 2.2mm



Cartridge Syringe  
DE-228 1.8mm  
DE-229 2.2mm



Euro Breech loading  
DE-230 1.8ml  
DE-231 2.2ml



EuroBreech loading  
DE-232 1.8ml  
DE-233 2.2ml



Manual Aspirating Syringe  
DE-234 1.8ml  
DE-235 2.2ml



Quatro  
DE-236 1.8ml  
DE-237 2.2ml




Trio  
DE-238 1.8ml  
DE-239 2.2ml



Swedish Type  
DE-240 1.8ml  
DE-241 2.2ml



Ject Ultima  
DE-242 1.8ml  
DE-243 2.2ml  
Injects 0.05 ml with every click Supplied with  
angled and straight tubes



Ject angled  
DE-244 1.8ml  
DE-245 2.2ml



Ject straight  
DE-246 1.8ml  
DE-247 2.2ml



Disc System  
DE-248 1.8ml  
DE-249 2.2ml



Disc System  
DE-250 1.8ml


# EXTRACTING FORCEPS SET

Извлечение щипцов  
Ensemble de pinces d'extraction  
• Zange sammeln

• Витягивающий комплект щипцов  
• Extracción del conjunto de fórceps  
• استخراج ملقط مجموعة



**DE-251**  
Extracting Forceps  
Set of 10 Pieces  
Extracting Forceps Fig. 2  
Extracting Forceps Fig. 7  
Extracting Forceps Fig. 13  
Extracting Forceps Fig. 17  
Extracting Forceps Fig. 18  
Extracting Forceps Fig. 22  
Extracting Forceps Fig. 33  
Extracting Forceps Fig. 51A  
Extracting Forceps Fig. 67A  
Extracting Forceps Fig. 79



**DE-252**  
Extracting Forceps  
Set of 10 Pieces  
Extracting Forceps Fig. 2  
Extracting Forceps Fig. 7  
Extracting Forceps Fig. 13  
Extracting Forceps Fig. 17  
Extracting Forceps Fig. 18  
Extracting Forceps Fig. 22  
Extracting Forceps Fig. 33  
Extracting Forceps Fig. 51A  
Extracting Forceps Fig. 67A  
Extracting Forceps Fig. 79



**DE-253**  
Extracting Forceps  
Set of 13 Pieces  
Extracting Forceps Fig. 2  
Extracting Forceps Fig. 7  
Extracting Forceps Fig. 13  
Extracting Forceps Fig. 17  
Extracting Forceps Fig. 18  
Extracting Forceps Fig. 22  
Extracting Forceps Fig. 33  
Extracting Forceps Fig. 51A  
Extracting Forceps Fig. 67A  
Extracting Forceps Fig. 79  
Extracting Forceps Fig. 86C  
Extracting Forceps Fig. 89  
Extracting Forceps Fig. 90



**DE-254**  
Extracting Forceps Children Pattern  
Set of 10 Pieces  
Extracting Forceps Fig. 13S  
Extracting Forceps Fig. 22S  
Extracting Forceps Fig. 30S  
Extracting Forceps Fig. 33S  
Extracting Forceps Fig. 37  
Extracting Forceps Fig. 38  
Extracting Forceps Fig. 39  
Extracting Forceps Fig. 39L  
Extracting Forceps Fig. 39R  
Extracting Forceps Fig. 51S



**DE-255**  
Extracting Forceps English adult Adult  
Set of 10 Pcs



**DE-256**  
Extracting Forceps American Adult  
Set of 10 Pcs




**DE-257**  
Extracting Forceps American Adult  
Set of 8 Pcs



**DE-258**  
Extracting Forceps Pedro  
Set of 7 Pcs



**Thumb Handle**



**Shortblast Handle**



**Half Gold Handle**



**Round With Edge Handle**



**Thom with Flat Handle**



**Finger Thom Handle**



**Anatomical Handle**



**Eye Style Handle**





# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

## Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.


### Apexo Root Elevators




### Luxator Root Elevators



### Bein Root Elevators



### Seldin Root Elevators



## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

**FUNCTION:** To loosen tooth or root from bony socket prior to placement of the extraction forceps.

### Laxutor Root Elevators

Long Type



Laxutor  
Root Elevator



**DE-318**  
Laxutor  
1.5mm  
Str.



**DE-319**  
Laxutor  
2mm  
Str.



**DE-320**  
Laxutor  
3mm  
Str.



**DE-321**  
Laxutor  
5mm  
Str.



**DE-322**  
Laxutor  
3mm  
Cvd.



**DE-323**  
Laxutor  
5mm  
Cvd.



**DE-324**  
Laxutor  
3mm  
Upper Angle



**DE-325**  
Laxutor  
3mm  
Double Angle



**DE-326**  
Laxutor  
3mm  
Str.



**DE-327**  
Laxutor  
5mm  
Str.




# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

## Laxutor Root Elevators



lower roots elevator



Upper root elevator with helicoidal tip


## Luxator Roots Elevator



### Function of Root Elevators:


To loosen tooth or root from bony socket prior to placement of the extraction forceps.

## Apical Root Elevators



Apical

## Apical Root Elevators



Apical




## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

**FUNCTION:** To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

### Cryer Root Elevators



### Seldin Root Elevators




# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber


- Слеватор корневой экстракционный
- Elevadores para raíces
- المصاعد الجذر

**FUNCTION:** To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

## Apical Root Elevators




## Apical Root Elevators



## Flohr Root Elevators




## Heidbrink Root Elevators



## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий экстракційний
- Elevadores para raíces
- المصاعد الجذر



**FUNCTION:** To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

### ENGLISH PATTERNS






# ROOT ELEVATOR


- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

## Lindo Levian Root Elevators



### Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

### Friedmann Roots Elevator



**DE-402**  
Friedmann 4F  
Str.

**DE-403**  
Friedmann 5F  
Left

**DE-404**  
Friedmann 6F  
Right

### Luxating Elevator With Periotome Tip



**DE-405**  
3mm Luxating  
Elevator / Periotome

**DE-406**  
4mm Luxating  
Elevator / Periotome

**DE-407**  
5mm Luxating  
Elevator / Periotome

### Wisdom Roots Elevator



**DE-408**  
Wisdom  
for 18 and 38

**DE-409**  
Wisdom  
for 28 and 48

**Function of Root Elevators:** To loosen tooth or root from bony socket prior to placement of the extraction forceps.

### Bayonet Roots Elevator Set



**DE-410**  
Left

**DE-411**  
Right

**DE-412**  
bayonet

**DE-413**  
Curved


**DE-414**  
Straight

# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber


- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

## Miller Root Elevators



**Function of Root Elevators:** To loosen tooth or root from bony socket prior to placement of the extraction forceps.

## Claw-Lift Roots Elevator



## Schlemmer Roots Elevator





## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

**FUNCTION:** To Loosen Tooth or Root From Bony Socket Prior to Placement of The Extraction Forceps.

**Laxutor Root Elevators Set of 9 Pcs**



DE-424  
2mm Str.



DE-425  
3mm Str.



DE-426  
5mm Str.



DE-427  
3mm Cvd.



DE-428  
5mm Cvd.



DE-429  
3mm Angle




DE-430  
3mm Angle



DE-431  
3mm Str.



DE-432  
5mm Str.



# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber


- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

## Hylin Roots Elevator



**Function of Root Elevators:** To loosen tooth or root from bony socket prior to placement of the extraction forceps.

## Kopp Roots Elevator



## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

### Luxier Roots Elevator



DE-447  
Luxier  
Cvd.

DE-448  
Luxier  
Str.



### Heber Roots Elevator



DE-449  
Heber  
Bayonet

DE-450  
Heber  
Bayonet



### Bein Roots Elevator



DE-451  
Bein,  
3mm

DE-452  
Bein,  
4mm

**Function of Root Elevators:** To loosen tooth or root from bony socket prior to placement of the extraction forceps.

### Bein Roots Elevator



DE-454  
2mm  
straight  
serrated

DE-455  
3mm  
straight

DE-456  
4mm  
Curved

DE-458  
2mm  
straight  
serrated

DE-459  
3mm  
straight

DE-460  
4mm  
Curved

DE-462  
2mm  
straight  
serrated

DE-463  
3mm  
straight

DE-464  
4mm  
Curved

DE-453  
Bein  
2mm, curved  
serrated

DE-457  
Bein  
2mm, curved  
serrated

DE-461  
Bein  
2mm, curved  
serrated



# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

## W. Barry Roots Elevator



**Function of Root Elevators:** To loosen tooth or root from bony socket prior to placement of the extraction forceps.



## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر



DE-482  
Fig. 1  
Chompret



DE-483  
Fig. 2  
Chompret



DE-484  
Fig. 3  
Chompret



DE-485  
Fig. 4  
Chompret



*Lecluse elevator for lower third molar dislocation*



*Chompret syndesmotome, sickle pattern*



DE-486



DE-487



DE-488



DE-489



DE-490



DE-491



DE-492

### Root Pick Elevator With Gold Flexible Points

# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

## Heidbrink Roots Elevator



**DE-493**  
Heidbrink  
Straight

**DE-494**  
Heidbrink  
Left

**DE-495**  
Heidbrink  
Right

**DE-496**  
Heidbrink  
Straight

**DE-497**  
Heidbrink  
Left

**DE-498**  
Heidbrink  
Right

## Syndesmotomes Roots Elevator



**DE-499**  
Straight

**DE-500**  
Left

**DE-501**  
Right

## Periotome Roots Elevator

## Syndesmotomes



**DE-502**  
Periotome  
2.1mm  
Straight

**DE-503**  
Periotome  
2.6mm  
Straight

**DE-504**  
Periotome  
2.1mm  
Angular

**DE-505**  
Periotome  
2mm  
Curved

**DE-506**  
Periotome  
2.0mm

**DE-507**  
Periotome  
2.5mm

**DE-508**  
Periotome  
3.0mm

**DE-509**  
Periotome  
4.0mm




## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

**FUNCTION:** To loosen small root fragments from bony socket

### Heidbrink Roots Elevator

### Warwick-James



## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر



**DE-519**  
1S  
(Straight 1mm)



**DE-520**  
2S  
(Straight 2mm)



**DE-521**  
3S  
(Straight 3mm)



**DE-522**  
5S  
(Straight 5mm)



**DE-523**  
3C  
(Curved 3mm)



**DE-524**  
5C  
(Curved 5mm)



**DE-525**  
3CA  
(Curved Angle 3mm)



**DE-526**  
3IC  
(Inverted Curved 3mm)



**DE-527**  
DE 3  
(Dual Edge 3mm / 1.5mm)




**DE-528**  
DE 5  
(Dual Edge 5mm / 3mm)

## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر



**DE-529**  
Luxators Set of 7 Pcs



**DE-530**  
Luxators Set of 7 Pcs



**DE-531**  
Luxators Set of 8 Pcs



**DE-532**  
Root Elevators Set of 10 Pcs



**DE-533**  
Elevators Set of 8 Pcs



**DE-534**  
Elevators Set of 8 Pcs




# ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Слеватор корневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

## PDL ELEVATORS



## ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

## PDL ELEVATORS

Delicate Arrow tip



4,0mm straight  
4,0mm curved

Luxivators straight  
atraumatic luxation of the tooth



2,5mm straight  
3,5mm straight  
4,5mm straight

Luxivators curved  
atraumatic luxation of the tooth




2,5mm curved  
3,5mm curved  
4,5mm curved



DE-537  
PDL Elevators  
Set of 8 Pcs

## Our Scalers

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 18 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. of course, any manufacturer can make these claims, but we guarantee and certify ours.



- High-polish mirror finish tips.
- Manufactured within the guidelines of the FDA quality assurance system for good manufacturing practices.
- Highest grade German stainless steel, razor sharp, durable tips. Rockwell hardened to 60 RHc.
- 27 of the most popular designs from graceys to McCalls to columbias, and many european patterns.
- All instruments are guaranteed against defects in material and workmanship.

- All instruments bear the ce mark, emblematic of the european community's highest standards.
- Customer satisfaction guarantee, too!
- Less than 20 grams ... Lightweight for less hand fatigue.
- Manufactured to the most exacting standards so the 50th instrument you buy is exactly like your first.
- You decide ... Why pay more for the same high-quality instruments?



## VALUED BUYERS

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours. We regularly manufacturing below articles in showing handles & pattern A to L for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in materials, sizes, shapes and patterns.



A-Hollow Handle Ø 10 mm



B-Hollow Handle Ø 10 mm



C-Hollow Handle Ø 10 mm



D-Hollow Handle Ø 10 mm



E-Hollow Handle Ø 10 mm



F-Hollow Handle Ø 10 mm



G-Hollow Handle Ø 10 mm



H-Hollow Handle Ø 6.4 mm



I-Hollow Handle Ø 8 mm



J-Hollow Handle Ø 10 mm



K-Soild Handle Ø 6 mm



L-Soild Handle Ø 6 mm

# MOUTH MIRRORS

- зеркал
- Espejos de la boca
- Mundspiegel
- дзеркал
- Miroirs à la bouche
- مرايا الفم

## Plain



**DE-600**  
Fig. 3  
Ø20mm

**DE-601**  
Fig. 4  
Ø22mm

**DE-602**  
Fig. 5  
Ø24mm

**DE-603**  
Fig. 6  
Ø26mm



## Cone-Socket Plain



**DE-604**  
Fig. 3  
Ø20mm

**DE-605**  
Fig. 4  
Ø22mm

**DE-606**  
Fig. 5  
Ø24mm

## Magnifying



**DE-607**  
Fig. 3  
Ø20 mm

**DE-608**  
Fig. 4  
Ø22 mm

**DE-609**  
Fig. 5  
Ø24 mm



## Cone-Socket Magnifying



**DE-610**  
Fig. 3  
Ø20mm

**DE-611**  
Fig. 4  
Ø22mm

**DE-612**  
Fig. 5  
Ø24mm

## Rhodium



**DE-613**  
Fig. 3  
Ø20mm

**DE-614**  
Fig. 4  
Ø22mm

**DE-615**  
Fig. 5  
Ø24 mm



## Cone-Socket Rhodium



**DE-616**  
Fig. 3  
Ø20mm

**DE-617**  
Fig. 4  
Ø22mm

**DE-618**  
Fig. 5  
Ø24mm

## Front Surface



**DE-619**  
Fig. 3  
Ø20mm

**DE-620**  
Fig. 4  
Ø22mm

**DE-621**  
Fig. 5  
Ø24mm

**DE-622**  
Fig. 6  
Ø26mm



## Cone-Socket Front Surface



**DE-623**  
Fig. 3  
Ø20mm

**DE-624**  
Fig. 4  
Ø22mm

**DE-625**  
Fig. 5  
Ø24mm

## HANDLES FOR MOUTH MIRRORS

- Ручки для зеркал
- Manches pour miroirs à bouche
- Mundspiegelgriffe
- Ручки для дзеркал
- Mangos para espejos de boca
- مقابض لمرآيا الفم



**DE-626**  
Hollow Mirror  
Handle Round  
Ø 10 mm



**DE-627**  
Hollow Mirror  
Hexagonal  
Handle  
Ø 8 mm



**DE-628**  
Hollow  
Ø 8 mm



**DE-629**  
Hollow Mirror  
Hexagonal  
Handle  
Ø 10 mm



**DE-630**  
Hollow Mirror  
Handle  
Ø 8 mm



**DE-631**  
Hollow Mirror  
Handle  
Ø 10 mm



**DE-632**  
Mirror Handle  
Solid  
Ø 6 mm



**DE-633**  
Mirror Handle  
Solid  
Ø 6 mm



**DE-634**  
Mirror Handle  
Solid Octagonal  
Ø 6 mm



**DE-635**  
Mirror Handle  
Stem  
with Rule  
35mm



**DE-636**  
Mirror Handle  
Solid  
Stem with  
Rule 35mm



**DE-637**  
Cone-Socket  
Solid Octagonal  
Ø 6 mm



**DE-638**  
Cone-socket  
Solid  
Ø 6 mm



**DE-639**  
Cone-socket  
Hexagonal  
Handle  
Ø 6 mm

Cone-Socket Mirror Handle



# EXPLORERS

- Зонды
- Sondes
- Zahnsonden

- Зонди
- Sondas
- المستكشفين

Examination kit Set of 3 Pieces



DE-640  
Probe

DE-641  
Mirror

DE-642  
College Forcep  
15cm



Examination Set of 4 Pieces



DE-643  
Exp,23/17-A

DE-644  
Excavators

DE-645  
Tweezers

DE-646  
Mouth Mirror

DE-647

First Aid Examination Set of 3 Pieces



DE-648

First Aid Examination Set of 3 Pieces



DE-649

First Aid Examination Set of 3 Pieces



DE-650

First Aid Examination Set of 3 Pieces



## TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



**DE-651**  
Surgical  
11.5cm



**DE-652**  
London College  
15cm



**DE-653**  
Semken Taylor  
Str.



**DE-654**  
Semken Taylor  
Cvd.



**DE-655**  
Meriam  
16cm



**DE-656**  
Dressing  
14cm



**DE-657**  
Semken Taylor  
1x2 Straight



**DE-658**  
Semken Taylor  
1x2 Curved



**DE-659**  
Gilles  
15cm 1x2



**DE-660**  
14cm 1x2



**DE-661**  
Perry  
13cm



**DE-662**  
Atraumatic  
15cm Cvd



**DE-663**  
Atraumatic  
15cm Str.



**DE-664**  
College with Lock  
15cm Cvd.



**DE-665**  
Gerald  
17.5cm



**DE-666**  
with lock  
15cm



**DE-667**  
Atraumatic  
17.5cm

## TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



**DE-668**  
Adson-TC  
12cm



**DE-669**  
Adson-TC  
12cm 1x2



**DE-670**  
Adson-TC  
15cm 1x2



**DE-671**  
Adson  
12cm



**DE-672**  
Adson  
12cm 1x2



**DE-673**  
Adson  
15cm 1x2



**DE-674**  
Adson Brown  
12cm



**DE-675**  
Crane Kaplan  
15cm



**DE-676**  
Crane Kaplan  
15cm



**DE-677**  
Crane Kaplan  
15cm



**DE-678**  
Hemingway  
15cm



**DE-679**  
Miller  
15cm Str.



**DE-680**  
Miller  
15cm Cvd.



**DE-681**  
Dressing  
Tweezer  
14.5cm



**DE-682**  
College  
Curved  
Tip 15cm



**DE-683**  
Maryam  
Tweezers  
15cm



**DE-684**  
College  
Tweezer  
15cm



**DE-685**  
Perry  
Tweezer  
12.5cm



**DE-686**  
Self-locking  
Tweezer  
15cm



## EXPLORERS


- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

Single Ended Probe Solid Handle 4.5mm Length 150mm



### Explorers-Single-End


### Double-Ended Explorers



- Зонды
- Sondes
- Zahnsonden

- Зонди
- Sondas
- المستكشفين


Probes Double Ended 6mm, length 170mm



Pointed Probes Double Ended Handle 6mm, length 170mm



**DE-710**  
Williams



## EXPLORERS

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

EXPLORERS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

### PLUGGERS






# PERIODONTAL PROBES

- Зонды пародонтологические
- Sondes paradontales
- sonda parodontali
- зонди пародонтологічні
- Sondas paradontales
- المسابير اللثة

## SPREADERS

## HEAT CARRIERS




## PERIODONTAL PROBES

- Зонды пародонтологические
- зонди пародонтологічні
- Sondes parodontales
- Sondas parodontales
- sonde parodontali
- المسابر اللثة


### ENDODONTICS INSTRUMENTS

#### ANTERIOR



#### LUKS







## GRACEY MINI FIVE

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



**DE-777**  
Fig. 3/4 Gracey 3M



**DE-778**  
Fig. 5/6 Gracey 3M



**DE-779**  
Fig. 7/8 Gracey 3M  
Buccal and lingual/  
Palatal surfaces



**DE-780**  
Fig. 11/12 Gracey 3M  
Mesial



**DE-781**  
Fig. 13/14 Gracey 3M  
Distal

**DE-776**  
Fig. 1/2  
Gracey 3m



## GRACEY AFTER FIVE

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



**DE-783**  
Gracey 3/4 AF



**DE-784**  
Gracey 5/6 AF



**DE-785**  
Gracey 7/8 AF



**DE-786**  
Gracey 11/12 AF



**DE-787**  
Gracey 13/14 AF

**DE-782**  
Gracey  
1/2 AF



### Gracey Curettes



3...1	1...3
3...1	1...3

**DE-788** Fig. 1/2 GR - Gracey  
Buccal surfaces



5...1	1...5
5...1	1...5

**DE-789** Fig. 3/4 GR - Gracey  
Lingual/Palatal surfaces



5...1	1...5
5...1	1...5

**DE-790** Fig. 5/6 GR - Gracey



8...4	4...8
8...4	4...8

**DE-791** Fig. 7/8 GR - Gracey  
Buccal and lingual/  
palatal surfaces



8...4	4...8
8...4	4...8

**DE-792** Fig. 9/10 GR - Gracey  
Buccal and lingual/  
palatal surfaces more angled



8...4	4...8
8...4	4...8

**DE-793** Fig. 11/12 GR - Gracey  
Mesial



8...4	4...8
8...4	4...8

**DE-794** Fig. 13/14 GR - Gracey  
Distal






8...4	4...8
8...4	4...8

**DE-795** Fig. 15/16 GR - Gracey  
Mesial



8...4	4...8
8...4	4...8

**DE-796** Fig. 17/18 GR - Gracey  
Distal



## Sickle Scalers



DE-797 Fig. U15/33 - Towner-Jaquette



DE-798 Fig. U15/30 - Towner-Jaquette



DE-799 Fig. H5/33



DE-800 Fig. H6/H7



DE-801 Fig. 30/33 - Jaquette



DE-802 Fig. T2/T3 - Taylor



DE-803 Fig. 204




DE-804 Fig. 204S



DE-805 Fig. 11/12 - MC - McCall  
Mesial




DE-806 Fig. 13/14 - MC - McCall  
Distal




**Straight and curved sickle scalers**

Feature: 2 cutting edges curved facial component tipped ends



**Universal curettes**

are characterized by an angle of 90° between the terminal stem and the cutting component.



**Gracey curettes**

are characterized by an angle of 70° between the terminal stem and the cutting component.



### Universal Curettes



DE-807 Fig. 13S/14S - MC - McCall



DE-808 Fig. 17S/18S - MC - McCall  
Distal



DE-809 Fig. M23



DE-810 Fig. 204SD



DE-811 Fig. YG7/8-Younger Good



DE-812 Fig. L 1/2 - Langer



DE-813 Fig. L 3/4 - Langer




DE-814 Fig. L 5/6 - Langer



DE-815 Fig. M 23A




DE-816 Fig. M 23B




**Straight and curved sickle scalers**  
Feature: 2 cutting edges curved facial component tipped ends



**Universal curettes**  
are characterized by an angle of 90° between the terminal stem and the cutting component.



**Gracey curettes**  
are characterized by an angle of 70° between the terminal stem and the cutting component.



## Periodontal knives





0.9 mm



0.7 mm

Chisels  
Diamond coated margin trimmers



Thanks to the diamond-coated hand instruments the margins are evenly levelled.



## Margin trimmers



**DE-835**  
1TH2 Hirt



**DE-836**  
5TH6 Hirt



**DE-837**  
3TH4 Hirt



**DE-838**  
7TH8 Hirt




**DE-839**  
62NS61

2.0 mm



**DE-840**  
37NS36

1.5 mm




**DE-841**  
66NS65


1.0 mm



**DE-842**  
79CI180 Black



**DE-843**  
77CI178 Black



**DE-844**  
279CI280 Black



**DE-845**  
277CI278 Black



**DE-846**  
179CI180 Black



**DE-847**  
177CI178 Black

## CHISELS



**DE-848**  
1/2 Werdelstaedt



**DE-849**  
3/4 Werdelstaedt



**DE-850**  
7 Straight



**DE-851**  
5/6 Straight



**DE-852**  
7/10 Straight



**DE-853**  
8/9 Binangle



**DE-854**  
11/12 Binangle



**DE-855**  
40/41 Binangle


## Tartar Remover Instruments Solid Handle Ø 6 mm



**DE-856**  
Fig. 1



**DE-857**  
Fig. 1L  
Left



**DE-858**  
Fig. 1R  
Right



**DE-859**  
Fig. 2L  
left



**DE-860**  
Fig. 2R  
Right



**DE-861**



**DE-862**



**DE-863**




**DE-864**  
WHITE 28




**DE-865**  
WHITE 29



**DE-866**  
WHITE 31



**DE-867**  
WHITE 32



**DE-868**  
WHITE 33



**DE-869**  
WHITE 34



**DE-870**  
WESTCOTT 21



**DE-871**  
WESTCOTT 21B



**DE-872**



**DE-873**




**DE-874**  
WHITE



**DE-875**  
WHITE







## Ball Burnishers



## Hatches







Spatulas rounded and sharp




Spatulas Blade Thickness 0.22 mm



## Spatulas



## Spatulas-Pluggers





Water drop Our excavators have particularly sharp blades, free of burrs and perfectly rounded.



**DE-929**  
0.8 MM



**DE-930**  
0.8 MM



**DE-931**  
0.8 MM



**DE-932**  
1.1 MM



**DE-933**  
1.1 MM



**DE-934**  
1.2 MM



**DE-935**  
1.2 MM



**DE-936**  
1.4 MM



**DE-937**  
1.6 MM





**DE-938**  
1.9 MM


























DE-1005  
WILLIAMS



DE-1006  
WILLIAMS



DE-1007  
WILLIAMS



DE-1008



DE-1009



DE-1010



DE-1011




DE-1012



DE-1013



DE-1014







## Serrated Pluggers



Ø 0.9 MM		DE-1034 OP36M	Ø 1.2 MM
Ø 1.1 MM		DE-1035 OP36D	Ø 1.6 MM
Ø 1.6 MM		DE-1036 Op36	Ø 2.2 MM
Ø 1.0 MM		DE-1037 Op52	Ø 1.7 MM
Ø 2.2 MM		DE-1038 Zi51	Ø 2.6 MM
Ø 1.8 MM		DE-1039 Op44	Ø 2.2 MM
Ø 1.7 MM		DE-1040 1BR2	Ø 2.6 MM
Ø 2.6 MM		DE-1041 2BR3	Ø 3.6 MM





**DE-1042 MF2C**  
Fuzzi



**DE-1043**  
Pj4



**DE-1044 MB3C**  
Besek



**DE-1045**  
HZ3C



**DE-1046**  
CPRO3C



**DE-1047 Mo1**  
Monari



**DE-1048 TH28C**  
Hirt



**DE-1049**  
HZ25C



**DE-1050**  
Ck10



**DE-1051**  
CK28BC



**DE-1052 Op37**  
Carver



**DE-1053**  
28BR29

## Round excavators

### Cariology

- for the removal of softened, carious dentin
- for the removal of temporary fillings

### Surgery

- to sever the periodontal ligament prior to extraction
- for the removal of granulation tissue after extraction

Design: Round working end Curved cutting edge



Round Excavator, 1.0 mm  
DE-1054



Round Excavator, 1.5 mm  
DE-1055



Round Excavator, 2.0 mm  
DE-1056

## Triple-angled back action excavators



Triple-angled Back Action Excavator, 1.0 mm  
DE-1057



Triple-angled Back Action Excavator, 1.5 mm  
DE-1058



Triple-angled Back Action Excavator, 2.0 mm  
DE-1059



Triple-angled Back Action Excavator, 3.5 mm  
DE-1060



# FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)


- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLACK


PLUGGERS/CONDENSERS

Hollenback



MARQUETTE CONDENSER

MORTONSON



## FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

### RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM Back Action Filling Instruments

#### Tanner Filling Instruments




#### BURNISHERS

#### BEAVERTAIL

#### ACORN

#### LADMORE



DE-1073  
1T  
Tanner


## CARVERS INSTRUMENTS



- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente

- Инструменти для реставрації
- Instrumentos para modelar
- أدوات تقطيع

### P.K. THOMS WAXING INSTRUMENTS SOLID HANDLE 6mm




DE-1095  
Composite Instruments Blue  
Titanium with Cassettes



DE-1094  
P.k. Thomas Waxing Instruments Kit

### COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC





## CARVERS INSTRUMENTS

- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente
- Инструменти для реставрації
- Instrumentos para modelar
- أدوات تقطيع



DE-1102



DE-1103



DE-1104



DE-1105



DE-1106




Cassette For 5 Instruments  
DE-1107

## FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

### FILLING INSTRUMENTS (Gold Pointed) SOLID HANDLE 6MM, LENGTH 170MM



DE-1108 1.5mm / 1.5mm

DE-1109 2.0mm / 2.0mm

DE-1110 2.0mm / 2.0mm

DE-1111 2.5mm / 2.5mm

DE-1112 1.3mm / 0.9mm

DE-1113 1.6mm / 1.3mm

DE-1114 2.4mm / 1.8mm

DE-1115 2.4mm / 2.0mm

DE-1116 1.3mm / 1.5mm

DE-1117 1.8mm / 1.2mm

Filling instrument,  
spatula 6 mm solid handle,  
length 170 mm,

Filling instrument  
6 mm  
length 175 mm,


Filling instrument,  
ball 6 mm  
length 170 mm,

Filling instrument,  
ball / spatula  
6 mm  
length 175 mm,

Filling instrument,  
plan / plan 6 mm  
length 165 mm,

## COMPOSITE FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك



DE-1118

DE-1119

DE-1120

DE-1121

DE-1122

DE-1123

DE-1124

DE-1125

DE-1126



## COMPOSITE FILLING INSTRUMENTS



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

### Composite Gold Tip Instruments



## MISCELLANEOUS PRODUCTS

- Прочие товары
- Divers produits et services
- Sonstige Produkte

- Інші товари
- Varios Productos
- المنتجات المتنوعة



**DE-1138**  
X Ray Film Hangers



**DE-1139**  
X Ray Film Hangers



**DE-1140**  
X Ray Film Hangers



**DE-1141**  
X Ray Film Hangers



**DE-1142**  
X Ray Film Hangers



**DE-1143**  
X Ray Film Hangers



**DE-1144**  
X Ray Film Hangers



**DE-1145**



**DE-1146**



**DE-1147**



Dissecting Set  
**DE-1148** Set of 9 Pcs.  
**DE-1149** Set of 11 Pcs.



**DE-1150**  
Napkin Holder with Plastic Chain



**DE-1151**  
Napkin holder with metal chain



Mouth Props  
**DE-1152**  
Small

Mouth Props  
**DE-1153**  
Medium

Mouth Props  
**DE-1154**  
Large

### Key Rings



**DE-1155**



**DE-1156**



**DE-1157**



**DE-1158**



**DE-1159**




**DE-1160**

## CEMENT SPATULAS

- Шпатели для цемента
- Spatules à ciment
- Zementspatel

- Шпателі для цементу
- Espátulas para cemento
- ملاعق الأسمنت





## WAX KNIVES

- Ножи для воска
- Couteaux a cire
- Wachsmesser

- Ножі для воску
- Cuchillos para cera
- السكاكين الشمع

ALGINATE, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR



**DE-1183**  
Wooden Handle  
Fig. 5  
21.5 cm



**DE-1184**  
Wooden Handle  
Fig. 7  
19.5 cm curved



**DE-1185**  
Wooden Handle  
Fig. 8  
20 cm curved



**DE-1186**  
Wooden Handle  
Fig. 9  
22.5 cm curved



**DE-1187**  
Hook Knives



**DE-1188**  
Fig. 1  
Lessmann  
12.5 cm



**DE-1189**  
Fig. 2  
Lessmann  
17 cm



**DE-1190**  
Fig. 1  
with small spoon  
12.5 cm



**DE-1191**  
Fig. 2  
with large spoon  
17 cm



**DE-1192**  
Bacolite Handle  
Wax Knives,  
Small



**DE-1193**  
Bacolite Handle  
Wax Knives,  
Large



**DE-1194**  
Wooden Handle  
Wax Knives,  
Small



**DE-1195**  
Wooden Handle  
Wax Knives,  
Large



**DE-1196**  
Plaster Knife



**DE-1197**  
Fig. 1  
17 cm



**DE-1198**  
Fig. 3  
16 cm



**DE-1199**  
Fig. 4  
16 cm



**DE-1200**  
Special Inox  
19.5 cm

# ARTICULATORS

- Артикуляторы
- Artocilar
- Artikulations

- Артикулятори
- Articuler
- ارنیکو لاتورس



**DE-1201**

**Technical Datas:**

Intercondylar distance  
condylar inclination  
Bennet angle  
Incisal table inclination  
Distance between mounting plates  
Weight  
Opening (item 5000)  
Opening (item 5010)

120 mm.  
fixed at 30°  
fixed at 7°  
10°  
85 mm.  
400 gr.  
125°  
180°



**DE-1202**

**Technical Datas:**

Intercondylar distance  
condylar inclination  
Bennet angle  
Incisal table inclination  
Distance between mounting plates  
Weight  
Opening

110 mm.  
fixed at 30°  
fixed at 7°  
flat  
85 mm.  
450 gr.  
180°



**DE-1203**

Stainless Steel Finishing

**DE-1204**

Royal Blue



**DE-1205**



**DE-1206**



**DE-1207**



**DE-1208**



**DE-1209**



**DE-1210**



**DE-1211**





**DE-1212**




**DE-1213**



**DE-1214**  
Nakawalni  
Dental Laboratory  
Instrument







## RETRACTORS - MOUTH GAGS

- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكاشطات - الكمامات الفم



**DE-1251**  
HMG  
Heister  
130mm



**DE-1252**  
Mouth Gag Molt  
11 cm for (Children)



**DE-1253**  
Mouth Gag Molt  
14 cm for (Adults)



**Sternberg Retractors**  
**DE-1254** 140 mm  
**DE-1255** 180 mm



**Columbia Retractor**  
**DE-1256** 140 mm



**DE-1257**  
Wangen Holder  
vestibule  
15 cm, with or without  
lip protection,



**DE-1258**  
Wangen Holder  
vestibule 15 cm, with  
or without lip protection



**DE-1259**  
Wangen Holder  
with dual  
function, 17 cm



**DE-1260**  
140 mm

**DE-1261**  
150 mm

**Weider Retractor**

## RETRACTORS - MOUTH GAGS

- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكاشطات - الكمامات الفم



Weitlaner  
Fig. 1  
DE-1262 10.5cm  
DE-1263 13.5cm



Weitlaner  
Fig. 2  
DE-1264 10.5cm  
DE-1265 13.5cm



DE-1266  
Lip & Cheek Retractor  
Stainless Steel  
Double Ended



DE-1267  
Simplex Retractor  
Stainless Steel  
Single Ended



Weitlander  
DE-1268 11cm  
DE-1269 13cm  
DE-1270 15cm



Jennings  
DE-1271 11cm  
DE-1272 13cm  
DE-1273 15cm



DE-1274  
Denhard  
12 cm



DE-1275  
Roser-König  
16 cm



DE-1276  
Roser-König  
19 cm



## RETRACTORS

- Ранорасширители
- Ecarteurs
- Wundhaken
- Ранорозширювачі
- Separadores
- الكاشحات

### Wound, Bendable

Retractor, individually or as saving kit!



**DE-1277**  
Cheek Retainer  
One-sided  
Long



**DE-1278**  
Cheek Retainer  
One-sided  
Long



**DE-1279**  
Cheek Retainer  
One-sided  
Long



**Ergo Grip  
Retractor**  
**DE-1280** 24 cm  
ergonomic handle,



**Ivory Retractors**  
Width adjustable  
& fixable,  
**DE-1281** 20 cm



**Bruenings Retractors**  
tongue depressors  
**DE-1282** 19 cm



**Minnesota  
Retractors**  
**DE-1283** 16 cm



**Ribbon Retractors**  
**DE-1284** 20 cm  
easy to bend  
by hand,



**Bishop Retractors**  
retractor and tongue  
pushers Bishop with  
comfortable margin  
**DE-1285** 16 cm



**Combination  
Retractors**  
**DE-1286** 25 cm

## RETRACTORS

- Ранорасширители
- Ecarteurs
- Wundhaken

- Ранорозширювачі
- Separadores
- الكاشطات



**DE-1287**  
Unique Retractors  
Length 175 mm  
angled 45 °



**DE-1288**  
Unique Retractors  
Length 175 mm  
angled 90°



**DE-1289**  
Retractors  
16cm



**DE-1290**  
Retractors  
19cm



**DE-1291**  
Mini-Langenbeck  
Retractors  
16 cm, hooks 22 mm



**DE-1292**  
US army Retractors  
21 cm, 2



**DE-1293**  
ALM Retractor  
about 8 cm



**DE-1294**  
Senn Mueller Retractor  
17 cm, hooks 22 mm

## RETRACTORS

- Ранорасширители
- Ecarteurs
- Wundhaken
- Ранорозширювачі
- Separadores
- الكاشحات



**DE-1295**  
Langenbeck Retractors  
surgical-flat  
21 cm, hook  
30x16 mm



**DE-1296**  
Langenbeck Retractors  
surgical-flat  
21 cm, hook  
42x11 mm



**DE-1297**  
Langenbeck Retractors  
surgical-flat  
21 cm, hook  
63x12 mm



Langenbeck Retractors  
surgical-flat  
21 cm, hook  
**DE-1298** 30x16 mm  
**DE-1299** 42x11 mm  
**DE-1300** 63x12 mm



Langenbeck Retractors  
21 cm, hook  
**DE-1301** 30x11 mm  
**DE-1302** 30x16 mm



**DE-1303**  
Middelborpf  
Retractors  
21 cm



**DE-1304**  
Farabeuf  
(120 X 12mm  
Set of 2 pices)



**DE-1305**  
Retractor Seldin



**DE-1306**  
Feldman Retractor  
14cm



## GUM SCISSORS

- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren

- Ножиці для ясен
- Tijeras para encias
- مقص صمغ



**DE-1307**  
Iris  
115mm  
Straight



**DE-1308**  
Iris  
115mm  
Curved



**DE-1309**  
Goldman-Fox  
130mm  
Saw edge, straight



**DE-1310**  
Goldman-Fox  
130mm  
Saw edge, Curved



**DE-1311**  
La Grange  
115mm  
Saw edge, double Curved



Mayo  
**DE-1312** 140mm  
**DE-1313** 170mm  
**DE-1314** 230mm  
Straight



Mayo  
**DE-1315** 140mm  
**DE-1316** 170mm  
**DE-1317** 230mm  
Curved



Metzenbaum  
**DE-1318** 145mm  
**DE-1319** 180mm  
**DE-1320** 200mm  
Straight



Metzenbaum  
**DE-1321** 145mm  
**DE-1322** 180mm  
**DE-1323** 200mm  
Curved



**DE-1324**  
Kelly  
160mm  
Straight



**DE-1325**  
Kelly  
160mm  
Curved



**DE-1326**  
Iris  
115mm  
Straight



**DE-1327**  
Iris  
115mm  
Curved



**DE-1328**  
La Grange  
115mm  
Double Curved



**DE-1329**  
Goldman-Fox  
130mm, Straight



**DE-1330**  
Goldman-Fox  
130mm, Curved



**DE-1331**  
Kelly  
160mm, Curved



**DE-1332**  
Iris  
115mm, Straight



**DE-1333**  
Iris  
115mm, Curved



**DE-1334**  
Micro Iris  
90mm, Curved

## GUM SCISSORS

- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren

- Ножиці для ясен
- Tijeras para encias
- مقص صمغ



**DE-1335**  
Iris  
115mm  
Angled



**DE-1336**  
Iris  
115mm  
Curved



**DE-1337**  
Goldman-Fox  
130mm  
Saw edge,  
straight



**DE-1338**  
Goldman-Fox  
130mm  
Saw edge,  
curved




**DE-1339**  
Goldman-Fox  
130mm  
Saw edge,  
double curved



**DE-1340**  
La Grange  
115mm  
Saw edge,  
double curved



**DE-1341**  
Neumann  
120mm  
Saw edge,  
double curved



**DE-1342**  
Kilner (Ragnell)  
120mm  
Straight



**DE-1343**  
Kilner (Ragnell)  
120mm  
Curved



**DE-1344**  
Joseph  
140mm  
Straight



**DE-1345**  
Joseph  
140mm  
Curved



**DE-1346**  
Kelly  
160mm  
Straight



**DE-1347**  
Kelly  
160mm  
Curved



**DE-1348**  
Locklin  
160mm  
Saw edge,  
angular



**DE-1349**  
Locklin  
120mm  
Saw edge,  
angular



**DE-1350**  
Dean  
170mm  
Saw edge,  
angular

Saw edge, angular



**DE-1352**  
Spencer  
90mm



**DE-1353**  
Spencer  
130mm



**DE-1354**  
Spencer  
115mm



**DE-1355**  
Noyes  
110mm, Straight



**DE-1356**  
Noyes  
110mm  
Curved



**DE-1357**  
Noyes  
110mm  
Angular



**DE-1358**  
Castroviejo  
150mm  
Curved

## GUM SCISSORS

## NEEDLE HOLDERS

- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren

- Ножиці для ясен
- Tijeras para encias
- مقص صمغ

- Иглодержатели
- Porte-aiguilles
- Nadelhalter

- Голкотримачі
- Porta-agujas
- اصحاب إبرة



**DE-1359**  
130mm  
Straight



**DE-1360** Mayo 130mm  
**DE-1361** Mayo 160mm  
**DE-1362** Mayo 180mm  
Straight Curved



**DE-1363** Mayo 140mm  
**DE-1364** Mayo 200mm  
Straight Curved



**DE-1365** Metzenbaum 145mm  
**DE-1366** Metzenbaum 200mm  
Straight Curved



**DE-1367**  
Beebe 110mm  
Saw edge, straight




**DE-1368**  
Beebe 110mm  
Saw edge, curved



**DE-1369**  
Beebe 110mm  
Saw edge, straight



**DE-1370**  
Beebe 110mm  
Saw edge, curved



**DE-1371**  
Universal 120mm  
Saw edge



**DE-1372**  
Universal 120mm



**DE-1373** Ryder 130mm  
**DE-1374** Ryder 150mm  
**DE-1375** Ryder 180mm



**DE-1376**  
Halsey 130mm



**DE-1377**  
Baumgartner 130mm



**DE-1378**  
Derf 120mm



**DE-1379**  
Delicate 150mm



**DE-1380** Mayo-Hegar 130mm  
**DE-1381** Mayo-Hegar 160mm  
**DE-1382** Mayo-Hegar 180mm  
**DE-1383** Mayo-Hegar 200mm



**DE-1384**  
Crile-Wood 150mm  
straight




**DE-1385**  
Crile-Wood 150mm  
curved



## NEEDLE HOLDERS

- Иглодержатели
- Porte-aiguilles
- Nadelhalter

- Голкотримачи
- Porta-agujas
- اصحاب إبرة



Ryder  
DE-1386 130mm  
DE-1387 150mm  
DE-1388 180mm



DE-1389  
Halsey  
130mm



DE-1390  
Baumgartner  
130mm



DE-1391  
Derf  
120mm



DE-1392  
Delicate  
150mm



DE-1393  
Adson  
170mm



Mayo-Hegar  
DE-1394 130mm  
DE-1395 160mm  
DE-1396 180mm  
DE-1397 200mm



DE-1398  
Crile-Wood  
150mm



Castroviejo  
DE-1399 str 140mm cvd  
DE-1400 str 180mm cvd



Olsen-Hegar Needle  
Holder With Scissors  
DE-1401 140mm



Gillies Needle  
Holder With Scissors  
DE-1402 160mm



Mathieu  
DE-1403 140mm  
DE-1404 170mm




Castroviejo  
DE-1405 str 140mm cvd  
DE-1406 str 180mm cvd



Olsen-Hegar Needle  
Holder With Scissors  
DE-1407 140mm



Gillies Needle  
Holder With Scissors  
DE-1408 160mm



Hartmann (Baby  
Mosquito)  
DE-1409 90mm straight



Halstead Mosquito  
DE-1411 125mm  
Straight

DE-1410  
125mm  
Curved

DE-1412  
Curved

## GUM SCISSORS

- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren

- Ножиці для ясен
- Tijeras para encias
- مقص صمغ



**DE-1413**  
Goldman Fox Scissors  
12.5cm Curved



**DE-1414**  
Iris Scissors  
11.5cm Curved  
Serrated For Slip-Free Cutting



**DE-1415**  
LaGrange Scissors  
11.5cm Double Curved  
Serrated For Slip-Free Cutting

Double  
Curved



**DE-1416**  
Castroviejo Micro Scissors  
12cm Curved



**DE-1417**  
Super Cut Goldman Fox  
13.2cm Curved  
Serrated For Slip-Free Cutting



**DE-1418**  
Super Cut Iris  
11cm Straight  
Serrated For Slip-Free Cutting



**DE-1419**  
Super Cut LaGrange  
11cm Curved  
Serrated For Slip-Free Cutting



**DE-1420**  
Super Cut LaGrange  
11.5cm  
Double Curved  
Serrated For Slip-Free Cutting

Double  
Curved

## NEEDLE HOLDERS

- Иглодержатели
- Porte-aiguilles
- Nadelhalter

- Голкотримачи
- Porta-agujas
- اصحاب إبرة



**DE-1421**  
Lichtenberg-  
Ryder T.C  
14 cm



**DE-1422**  
Lichtenberg-  
Ryder T.C  
14 cm



**DE-1423**  
Lichtenberg-  
Ryder T.C  
17 cm



**DE-1424**  
Lichtenberg-  
Ryder T.C  
17 cm



**DE-1425**  
Barraquer  
18 cm



**DE-1426**  
Barraquer  
18 cm



**DE-1427**  
Barraquer T.C  
18 cm



**DE-1428**  
Barraquer T.C  
18 cm



**DE-1429**  
Tönnis  
18 cm



**DE-1430**  
Castroviejo  
13 cm  
Straight



**DE-1431**  
Kalt Needle Holders T.C  
14 cm Straight



**DE-1432**  
Peet Splinter Forceps  
11.5 cm

**DE-1433**  
La Grange  
Curved




**DE-1434**  
All-Purpose Pliers



# HAEMOSTATIC FORCEPS & TOWEL CLAMPS

- Зажимы артериальные, зажимы для салфеток
- Pincas hemostaticas, Pincas a compresas
- Arterienklemmen, Tuscheklemmen
- Затискачі артеріальні, затискачі для серветок
- Pinzas hemostaticas, Pincas para compresas
- المنقط هيموستاتيك المشابك منشقة




**DE-1436**  
curved

Halstead Mosquito 1x2  
**DE-1435** 125mm  
straight



**DE-1438**  
curved

Kelly  
**DE-1437** 140mm  
straight




**DE-1443** 140mm  
**DE-1444** 160mm  
**DE-1445** 180mm  
**DE-1446** 200mm  
curved

Rochester-Pean  
**DE-1439** 140mm  
**DE-1440** 160mm  
**DE-1441** 180mm  
**DE-1442** 200mm  
straight



**DE-1451** 140mm  
**DE-1452** 160mm  
**DE-1453** 180mm  
**DE-1454** 200mm  
curved

Rochester-Kocher 1x2  
**DE-1447** 140mm  
**DE-1448** 160mm  
**DE-1449** 180mm  
**DE-1450** 200mm  
straight



**DE-1456**  
curved

**DE-1455**  
Kelly-Rankin  
160mm  
straight



Lorna  
**DE-1457** 100mm  
**DE-1458** 130mm



**DE-1459**  
Allis-Baby 3:4  
150mm



**DE-1460**  
Allis 4:5  
200mm



**DE-1461**  
Allis 6:7  
200mm



Ball Socket  
**DE-1462** 140mm  
Not-perforating and  
not-damaging



**DE-1463**  
Backhaus  
130mm



**DE-1464**  
Backhaus  
85mm



**DE-1465**  
Jones  
90mm



**DE-1139**  
curved

Gross Maier  
Dressing Forceps  
**DE-1466** 200mm straight



Foerster  
Dressing Forceps  
**DE-1467** 200mm



Young  
Tongue Forceps  
**DE-1468** 170mm



Collins  
Tongue Forceps  
**DE-1469** 170mm

## RUBBER DAM CLAMPS FORCEPS

- Инструменты для системы коффердам
- Pincas porte-crampons a dique
- Kofferdam Klammerzangen
- Инструменти для системи коффердам
- Pinzas de engrapar para dique de caucho
- السد العطار الملقط المشابك



### DE-1470

On request the Content of the rubber dam kit can be changed with other instruments.

#### Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Denture Plate
- 05 Rubber Dam Frame
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case
- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7



**DE-1474**  
Brewer



**DE-1475**  
Stocks



**DE-1476**  
Palmer




**DE-1477**  
Ivory light weight



**DE-1478**  
Anisworth



**DE-1479**  
Anisworth



**DE-1480**  
Anisworth



**DE-1481**  
Ivory



**DE-1482**  
Seaparing  
Elastic Plier

## RUBBER DAM CLAMPS KITS

- Резиновая дамба зажимы комплекты
- Trousses de caoutchouc
- Gummidammklammern

- Гумова дамба затискачі комплекти
- Conjuntos de abrazaderas de retención de caucho
- المطاط السد المشابك مجموعات



**DE-1483**

Stainless Steel Gold Plated Rubber Dam Clamp Tray with 11 Different Clamps  
Fig. 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



**DE-1484**

Stainless Steel Rubber Dam Clamp Tray with 12 Different Clamps  
Fig. 209, 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7




**DE-1485**

Stainless Steel Rubber Dam Clamp Tray with 8 Different Clamps  
Fig. 1, 26, 29, 51, 60, 202, 205, 210



**DE-1486**

Stainless Steel Rubber Dam Clamp Tray  
(without clamps)




**DE-1487**

On request the Content of the rubber dam kit can be changed with other instruments.

**Content:**

- |                              |                              |
|------------------------------|------------------------------|
| 01 Ivory                     | 08 Rubber Dam Clamp Fig. 209 |
| 02 Ivory light weight        | 09 Rubber Dam Clamp Fig. 212 |
| 03 Clamps Stand              | 10 Rubber Dam Clamp Fig. 211 |
| 04 Rubber Dam Frame M        | 11 Rubber Dam Clamp Fig. 8   |
| 05 Rubber Dam Frame L        | 12 Rubber Dam Clamp Fig. 8A  |
| 06 Rubber Dam Box Lower Case | 13 Rubber Dam Clamp Fig. 14  |
| 07 Rubber Dam Box Upper Case | 14 Rubber Dam Clamp Fig. 14A |
|                              | 15 Rubber Dam Clamp Fig. 0   |
|                              | 16 Rubber Dam Clamp Fig. 00  |
|                              | 17 Rubber Dam Clamp Fig. 2   |
|                              | 18 Rubber Dam Clamp Fig. 2A  |
|                              | 19 Rubber Dam Clamp Fig. 7   |



**DE-1488**

Rubber Dum Template



## RUBBER DAM CLAMPS

- Кламмера для коффердама
- Crampons pour digue
- Kofferdam Klammern
- Кламмера для коффердама
- Grapas para dique de caucho
- المشابك السد المطاط



DE-1489  
Fig. 15



DE-1490  
Fig. 18



DE-1491  
Fig. 24



DE-1492  
Fig. 25



DE-1493  
Fig. 26



DE-1494  
Fig. 28



DE-1495  
Fig. 30



DE-1496  
Fig. 31



DE-1497  
Fig. 51



DE-1498  
Fig. 138



DE-1499  
Fig. 139



DE-1500  
Fig. 200



DE-1501  
Fig. 201



DE-1502  
Fig. 202



DE-1503  
Fig. 203



DE-1504  
Fig. 204



DE-1505  
Fig. 205



DE-1506  
Fig. 3



DE-1507  
Fig. W3



DE-1508  
Fig. 4



DE-1509  
Fig. 56A



DE-1510  
Fig. 7



DE-1511  
Fig. 7A



DE-1512  
Fig. 8



DE-1513  
Fig. 8A



DE-1514  
Fig. W8A



DE-1515  
Fig. 10



DE-1516  
Fig. 11



DE-1517  
Fig. 12A



DE-1518  
Fig. 13A



DE-1519  
Fig. 14



DE-1520  
Fig. 14A



DE-1521  
Fig. 14AN



DE-1522  
Fig. 22



DE-1523  
Fig. 27



DE-1524  
Fig. 207



DE-1525  
Fig. 29



DE-1526  
Fig. 206



DE-1527  
Fig. 208



DE-1528  
Fig. 209



DE-1529  
Fig. 1



DE-1530  
Fig. 1A



DE-1531  
Fig. 2



DE-1532  
Fig. 2A



DE-1533  
Fig. W2A



DE-1534  
Fig. 0



DE-1535  
Fig. 00



DE-1536  
Fig. 210



DE-1537  
Fig. 211



DE-1538  
Fig. 212



DE-1539  
Fig. 6



DE-1540  
Fig. 9



DE-1541  
Fig. N06



DE-1542  
Fig. U67



DE-1543  
Fig. 214



DE-1544  
Fig. P1



DE-1545  
Fig. P2



DE-1546  
Fig. G-1



DE-1547  
Fig. G-2

## CROWN REMOVER

- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen
- Инструменты для удаления коронок
- Arranca-corona
- ميزلات التاج



**DE-1548**  
Crown Spreaders



**DE-1550**  
Universal  
Fig. 0



**DE-1552** **DE-1553**



**DE-1554**

Manual Crown Removers  
Handle with 3 tips **DE-1551**  
Only 3 tips **DE-1552 to DE-1554**



**DE-1556** **DE-1557**



**DE-1558**

Manual Crown Removers  
Handle with 3 tips **DE-1555**  
Only 3 tips **DE-1556 to DE-1558**



**DE-1559**  
Manual Crown  
Remover



**DE-1560**  
Manual Crown Remover



**DE-1561**  
Complete with 3 tips  
and tightening keys



semi-automatic crown  
& bridge remover




**DE-1562**



**A**



**B**



**DE-1563**  
Manual Crown Remover



**A**



**B**



**C**



**D**



**E**

## MATRIX RETAINERS

- Матрицедержатели
- Porte Matrices
- Matrizenspanner

- Матрицеутримувачи
- Porta Matrices
- مصفوفة الخدم



DE-1564  
Fig. 5  
Ivory 8  
5 cm



DE-1565  
Fig. 6  
Ivory 8  
6 cm



DE-1566  
Fig. 7  
Ivory 8  
7 cm



DE-1567  
Fig. 1  
Tofflemire-Universal



DE-1568



DE-1569



DE-1570

### Siqveland Retainers



DE-1571  
Fig. 2  
Tofflemire-Junior



DE-1572  
Fig. 3  
Tofflemire-Senior



DE-1573  
Narrow



DE-1574  
Wide

### Amalgam Well



DE-1575  
Stainless Steel  
Amalgam Well



DE-1576  
Ivory 12



DE-1577  
Ivory 1



DE-1578  
Ivory




DE-1579  
Elliot



## AMALGAM INSTRUMENTS

- Инструменты для амальгамы
- Instruments pour amalgame
- Amalgamisinstrumente
- Инструменти для амальгами
- Instrumentos para amalgama
- الصكوك الملعفم



## AMALGAM INSTRUMENTS

- Инструменты для амальгамы
- Instruments pour amalgame
- Amalgamis\instrumente
- Инструменты для амальгами
- Instrumentos para amalgama
- الصكوك المنغم

### AMALGAM CARRIERS

#### TIP SIZES FOR DISTAL OR SINGLE-END

1.5 mm

2.0 mm

2.8 mm

1.5 mm

2.0 mm



**DE-1592**  
Mini/  
Single-End



**DE-1593**  
Regular  
Single-End



**DE-1594**  
Large  
Single-End



**DE-1595**  
Mini  
Distal



**DE-1596**  
Regular  
Distal

#### TIP SIZES FOR DOUBLE-END

1.5 mm

2.8 mm

3.2 mm



2.0 mm  
**DE-1597**  
Mini/  
Regular



2.0 mm  
**DE-1598**  
Large/  
Regular



2.0 mm  
**DE-1599**  
Jumbo/  
Regular



**DE-1600**  
Amalgam  
Well

**DE-1591**  
Single-End

#### TIP SIZES FOR DOUBLE-END

**DE-1601** - 1.5mm-2.0mm  
**DE-1602** - 2.0mm-2.8mm  
**DE-1603** - 2.0mm-3.2mm

**DE-1604** - 1.0mm-2mm  
**DE-1605** - 2.0mm-3mm  
**DE-1606** - 2.0mm-4mm

**DE-1607** - 1.5mm-2.0mm  
**DE-1608** - 2.0mm-2.8mm  
**DE-1609** - 2.0mm-3.2mm

**DE-1610** - 1.0mm-2mm  
**DE-1611** - 2.0mm-3mm  
**DE-1612** - 2.0mm-4mm

**DE-1613** - 1.5mm-2.0mm  
**DE-1614** - 2.0mm-2.8mm  
**DE-1615** - 2.0mm-3.2mm

**DE-1616** - 1.0mm-2mm  
**DE-1617** - 2.0mm-3mm  
**DE-1618** - 2.0mm-4mm

**DE-1619** - 1.5mm-2.0mm  
**DE-1620** - 2.0mm-2.8mm  
**DE-1621** - 2.0mm-3.2mm

**DE-1622** - 1.0mm-2mm  
**DE-1623** - 2.0mm-3mm  
**DE-1624** - 2.0mm-4mm

**DE-1625** - 1.5mm  
**DE-1626** - 2.0mm  
**DE-1627** - 2.0mm  
**DE-1628** - 3.0mm





## ROOT PLANNING FILES



**DE-1641**  
Krarner-Nevins  
Bone Chisel



**DE-1642**  
Krarner-Nevins  
Bone Chisel



**DE-1643**  
Molt 18.5cm  
Auto Clavable  
Plastic Tips.



**DE-1644**  
Mead 18.5cm  
Non Auto Clavable  
Plastic Tips.



**DE-1645**  
Miller



**DE-1646**  
Miller

## BONE FILES



**DE-1647**  
Miller 21



**DE-1648**  
Miller 1X  
cross cut



**DE-1649**  
Miller 2X  
cross cut



**DE-1650**  
Miller 52




**DE-1651**  
Wahl 1

# PERIODONTAL INSTRUMENTS

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie

- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة

## PERIOSTEEL ELEVATORS



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة

To separate tissue tooth or bone. to hold tissue away from surgical site.





## PERIODONTAL INSTRUMENTS



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die parodontologie
- Інструменти пародонтологічні
- Intrumentos para periodoncia
- الصكوك اللثة



DE-1689



DE-1690  
P24G



DE-1691  
titanium Periosteal Elevator for Implan5



DE-1692  
S1



DE-1693  
Hirschfeld



DE-1694  
9S



10mm

DE-1695  
Selden

13mm



4.8 mm

Prichard

10.6 mm

DE-1696  
Friedmann



6.9 mm

DE-1697  
Schmidlin

## SCALPEL HANDLES

- Ручки для скальпелей
- Manches de bistouris
- Skalpellgriffe
- Ручки для скальпелів
- Mangos de bisturi
- مشرط مقابض



**DE-1698**  
Fig. 3  
12.7cm



**DE-1699**  
Fig. 3G  
12.7cm



**DE-1700**  
Fig. 4  
13.7cm



**DE-1701**  
Double Handle  
No. 4 + No. 3  
16.3cm



**DE-1702**  
Fig. 7  
16.5cm



**DE-1703**  
Scalpel Handle Straight  
for blades #10-15  
14.5cm



**DE-1704**  
Micro Surgical Scalpel  
Handle titanium  
13.5cm



**DE-1705**  
Fig. 5  
14.5cm



**DE-1706**  
Fig. 5A  
14.5cm



**DE-1707**  
6 Way Adjustable Angled  
Scalpel Handle



**DE-1708**  
Metal tip  
220 mm  
Hammer With Wooden Handle



**DE-1709**  
Metal tip  
220 mm



**DE-1710**  
1  
For plaster



**DE-1711**  
2  
For metal



**DE-1712**  
Fret saw  
250 mm



**DE-1713**  
Small  
Saw Blades



**DE-1714**  
Large  
Saw Blades



**DE-1715**  
Ceri Saw



Special Inox  
**DE-1716** 170mm  
**DE-1717** 200mm



**DE-1718**  
Plaster Shears  
200 mm



Crucible Tongs  
**DE-1719** 250 mm  
**DE-1720** 320 mm



**DE-1721**  
Ceri Saw  
180 mm



**DE-1722**  
Blade  
0,12x5 mm  
Conf. 12 pcs



**DE-1723**  
Blade  
0,25x5 mm  
Conf. 12 pcs



**DE-1724**  
Blade  
0,25x3 mm  
Conf. 12 pcs





Behind Every Smile is a Story

- QUALITY - ЯКІСТЬ
- PRECISION - ТОЧНІСТЬ
- ACCURACY - ПРАВИЛЬНІСТЬ
- DURABILITY - ЗНОСОСТІЙКІСТЬ



# ORTHODONTIC INSTRUMENTS



- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**DE-1800**  
Height Gauge 0.18  
Aluminum



**DE-1801**  
Height Gauge 0.22  
Aluminum



**DE-1802**  
Height Gauge 0.18  
Stainless Steel



**DE-1803**  
Height Gauge 0.22  
Stainless Steel



**DE-1804**  
Moveable  
Height Gauge



**DE-1805**  
Endo Finger  
Ruler Span  
Measure Scale



**DE-1806**  
Endo Finger  
Ruler Span  
Measure Scale



**DE-1807**  
Endo Finger  
Ruler Span  
Measure Scale



**DE-1808**  
Positive  
Postioner



**DE-1809**  
Strip Saw



**DE-1810**  
Bone Bracket Gauge



**DE-1811**  
Endo Finger Ruler  
Span Measure Scale



**DE-1812**  
Endo Finger Ruler  
Span Measure Scale

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**DE-1813**  
**Measuring Loop**  
angled, length  
175 mm,  
measuring range  
0 - 40 mm



**DE-1814**  
**Measuring Loop**  
straight length  
175 mm,  
measuring range  
0 - 40 mm



**DE-1815**  
**Measuring loop, Castroviejo**  
angled, length 175 mm,  
measuring range 0 - 40 mm,  
with window



**DE-1816**  
**Measuring Loop**  
straight length 90 mm,  
measuring range  
0 - 40 mm,

**DE-1817**  
Castroviejo Calipers  
20mm  
90mm



**DE-1818**  
**Measuring Loop**  
angled, length 90 mm,  
measuring range  
0 - 40 mm,

**DE-1819**  
Castroviejo Calipers  
20mm, 45°  
90mm



**DE-1820**  
Castroviejo Calipers  
20mm  
90mm



**DE-1821**  
Fig. 1  
Iwanson  
spring caliper  
For metal and  
porcelain.



**DE-1822**  
Fig. 2  
Iwanson  
spring caliper  
For wax.



## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Інструменти ортодонтичні
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**DE-1823**  
Boley Gauge



**DE-1824**  
Vernier Gauge



**DE-1825**  
Aluminium Caliper  
up to 200mm



**DE-1826**  
Boley Gauge



**DE-1827**  
Delux Boley  
Gauge



**DE-1828**  
Beerendonk  
130mm  
Mini Gauge



**DE-1829**  
Compasso  
De Willis



**DE-1830**  
Dentaurum  
Ruler



**DE-1831**  
Ruler Wedge



**DE-1832**  
Ruler

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان

## Bracket Holder Tweezers



**DE-1833**  
Deluxe Bracket  
Holder Tweezers



**DE-1834**  
Deluxe Bracket  
Holder Tweezers



**DE-1835**  
Direct Bond Bracket  
Holder Tweezers



**DE-1836**  
Bracket Tweezers  
Slot Aligner



**DE-1837**  
Direct Bond  
Bracket Holder  
Tweezers, Small Tip



**DE-1838**  
Exact Bracket  
Holder Tweezers



**DE-1839**  
Bracket Placing  
Tweezers



**DE-1840**  
Easy Access  
Tweezers



**DE-1841**  
Bracket Placer  
Tweezers  
(Extra Fine Tip)



**DE-1842**  
Deluxe Bracket  
Holder Tweezers  
135mm



**DE-1843**  
Easy Eccess  
Tweezers



**DE-1844**  
Slipper  
Tweezers



**DE-1845**  
Direct Bonding  
Bracket Holders  
135mm



**DE-1846**  
Direct Bonding  
Bracket Holders  
136mm



**DE-1847**  
Direct Bonding  
Bracket Holders  
135mm



**DE-1848**  
Deluxe Bracket  
Holder



**DE-1849**  
Direct Bonding  
Bracket Holders  
140mm

# ORTHODONTIC INSTRUMENTS



- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente

- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**DE-1850**  
Band seater  
auto cleavable



**DE-1851**  
Band seater  
auto cleavable  
Triangle Shape



**DE-1852**  
Band seater  
auto cleavable  
Square Shape



**DE-1853**  
Band seater  
auto cleavable  
Round Shape



**DE-1854**  
Mershon Band  
Pusher  
140mm



**DE-1855**  
Band seating &  
removing  
instruments



**DE-1856**  
Ortho Light double  
ended  
ligature director  
scaler



**DE-1857**  
Ortho Light Band  
Pusher Scaler



**DE-1858**  
Band Seating &  
Removing  
Silicon Grip



**DE-1859**  
Ligature Instrument  
Silicon Grip



**DE-1860**  
Elastic Instrument  
Silicon Grip



**DE-1861**  
Distal Bender  
Silicon Grip



**DE-1862**  
Distal Bender  
Silicon Grip



**DE-1863**  
Distal Bender  
Silicon Grip



## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**DE-1864**  
Ortho light scaler double ended



**DE-1865**  
Ortho light scaler double ended



**DE-1866**  
Ortho light scaler double ended



**DE-1867**  
Twiel-on elastic module instrument



**DE-1868**  
Twiel-on elastic module instrument



**DE-1869**  
Schure Band Seater



**DE-1870**  
Ortho light scaler



**DE-1871**  
Ortho light scaler



**DE-1872**  
Ortho light scaler heavy



**DE-1873**  
Ortho instrument



**DE-1874**  
Ortho instrument

## ORTHODONTIC INSTRUMENTS




- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**DE-1875**  
Ortho Light double ended  
ligature director



**DE-1876**  
Ortho Light double ended ligature scaler director tivate coil spring  
are remove modules with one instruments



**DE-1877**  
Ortho band pusher scaler



**DE-1878**  
Ortho ligature director single ended



**DE-1879**  
Ortho Light double ended ligature  
director scaler



**DE-1880**  
Ortho Light double ended ligature scaler/director  
tivate coil spring are remove modules with one instruments



**DE-1881**  
Ortho Light double ended ligature director



**DE-1882**  
Ortholite Ligature Director Scaler



**DE-1883**  
Ortho light ligature director double ended



**DE-1884**  
Double ended cement spatula



**DE-1885**  
Ortho light scaler single ended

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان برومبشيس

### BOX JOINT PLIERS

#### DE-1886

Universal Dolphin Style  
Wire Bending & cutting pliers



#### DE-1892

Bracket Removing Plier curved  
Bracket Removing Plier curved for  
easy removal of brackets in the  
posterior region. Length 130mm



#### DE-1887 Band Removing Plier

The curved blade of the band  
removing plier can strip easily  
and gently remove directly to  
the tooth. Length 130mm



#### DE-1893

Young Wire Bending &  
Loop Forming Pliers



#### DE-1888

Tweed Loop  
forming Pliers



#### DE-1894

Wire Bending &  
Loop Forming Pliers



#### DE-1889

How Plier curved  
Universally usable tool for  
the multiband technology for  
wires up to .028".  
Length 130mm



#### DE-1895

Adams, small wire bending  
and loop forming Pliers



#### DE-1890 How Plier

Universally usable tool for the  
multiband technology for wires  
up to .028".  
Length 135mm



#### DE-1896

Bracket Removing Plier Round Box  
Bracket removing plier for comfortable  
and quick removal of the brackets.  
Length 135mm



#### DE-1891

Nance Loop  
forming Pliers



#### DE-1897

Weingart Plier medium  
Weingart plier that are ideal due to  
their curved tip for setting sheets and  
is capable of bending the ends.  
Length 145mm





## PLIERS FOR ORTHODONTICS & PROSTHETICS



- Щипцы для ортодонтии и протезирования
- Пинцес pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروتستیس

### BOX JOINT PLIERS

#### DE-1898 Adams Plier Round Box

Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



#### DE-1904 Adams Plier

Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



#### DE-1899

Falcon Beak Contouring Pliers.



#### DE-1905

Tweed Loop Formin Pliers



#### DE-1900

wires up to 0.020" / 0.7mm.



#### DE-1906 Bird Beak Plier Standard

Bird Beak Wire bending plier without notches, which is suitable as a universal utility tool. Length 125mm



#### DE-1901

Johnson Contouring Pliers



#### DE-1907

Mini Ligature cutters



#### DE-1902 Aderer

Large wire bending and loop forming pliers



#### DE-1908 Small

Fluch-cut distal and cutting Pliers with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



#### DE-1903

Tweed / angle rectangular arch forming & Torquing pliers



#### DE-1909 Large

Fluch-cut distal and cutting Pliers with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومبلیتیس

### BOX JOINT PLIERS

**DE-1910**

**Tweed / Angel**

1.78mm (.070") wide tips arch forming and torquing pliers



**DE-1916**

**Contouring Plier**



**DE-1911**

**Tweed / Angel**

1.78mm (.070") wide tips arch forming and torquing pliers



**DE-1917**

**Bands & Crowns Plier**  
in order to skirt



**DE-1912**

**Lingual / Palatal arch**  
forming pliers



**DE-1918**

**Ligature Forming Pliers**



**DE-1913**

**Aderer**

Facebow loop forming Pliers



**DE-1919**

for double back and triple back bands in 0.30" / 0.76mm & 0.36" / 0.91mm lingual arch wires



**DE-1914**

**Waldsavs wire bending and cutting pliers**



**DE-1920**

**Bird Beak Wire Bending & cutting pliers**



**DE-1915 Aderer Plier**

For the manufacturing of appliances and activation of wires for wires up to .030" .fine tip and up to .036" mean tip.  
Length 130mm



**DE-1921**

**Weingart Plier medium**

Weingart plier that are ideal due to their curved tip for setting sheets and is capable of bending the ends.  
Length 145mm





# PLIERS FOR ORTHODONTICS & PROSTHETICS



- Щипцы для ортодонтии и протезирования
- Пинсес pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- ڪمائيءَ لاءِ ٽٽوڙيل اسٽين بروسٽيٽيس

## L-KEY JOINT PLIERS



### DE-1922

Distal end cutter with safety hold, mini head, long handles.  
fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



### DE-1923

Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" to .025")



### DE-1924

Soft wire and ligature cutter, Mini fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



### DE-1925

Soft wire and ligature cutter, standard fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



### DE-1926

Soft wire and ligature cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



### DE-1927

Soft wire and ligature cutter, 45° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



### DE-1928

Hard wire cutter, straight Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### DE-1929

Hard wire cutter, 15° angled Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### DE-1930

Hard wire cutter, 15° angled Non-autoclavable

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومبلیتیس

## L-KEY JOINT PLIERS



**DE-1931**  
Band removing pliers



**DE-1932**  
Bracket removing pliers angled



**DE-1933**  
Bracket removing pliers straight



**DE-1934**  
Bracket removing pliers with plastic tip



**DE-1935**  
Adhesive removing pliers



**DE-1936**  
Weingart-mini  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



**DE-1937**  
Weingart  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



**DE-1938**  
How Straight  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



**DE-1939**  
Angled  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



**DE-1940**  
Wire bending and loop forming pliers  
Angled (Bird beak)

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")

# PLIERS FOR ORTHODONTICS & PROSTHETICS



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروتستيس

## L-KEY JOINT PLIERS



**DE-1941**  
Multi-purpose pliers, bird beak,  
De la rosa and cutter in one

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.56 mm (.018" x .022")



**DE-1946**  
Stop pliers V-Bebd - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.43 mm x 0.56 mm (.017" x .028")



**DE-1942**  
Cinch Back / tie back  
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**DE-1947**  
Lingual / palatal arch forming plies

- 0.76 mm (.030") / 0.91 mm (.036")



**DE-1943**  
Hammer head  
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**DE-1948**  
Torquing pliers set



**DE-1949**  
Tweed / angel  
Rectangular arch forming and  
torquing pliers

- 0.30 mm to 0.71 mm (.012" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**DE-1944**  
Nance  
Loop forming pliers

- 0.36 mm to 0.71 mm (.014" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**DE-1950**  
Surgical ball hook crimping pliers,  
angled



**DE-1945**  
Detailing step pliers / Bayonet  
bending pliers - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**DE-1951**  
Surgical ball hook crimping  
pliers, straight



- Щипцы для ортодонтии и протезирования
- Pincettes pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان وبروتيسيس

## L-KEY JOINT PLIERS




**DE-1952**  
Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/.76mm.



**DE-1958**  
Falcon beak contouring plier.




**DE-1953**  
Bird beak plier for bending round, square and rectangle wire up to .30"/.76mm.



**DE-1959**  
For placement and removal of arch wires, pins and other auxiliaries.




**DE-1954**  
Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.



**DE-1960**  
For placement and removal of arch wires, pins and other auxiliaries in angle design



**DE-1955**  
Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.




**DE-1961**  
contouring Plier.



**DE-1962**  
Bending plier. In order to skirt molars bands. With grooved. Recommended wire. 016", 018" and 022" 0.41mm, 0.46mm and 0.56mm



**DE-1956**  
Adams pliers for hooks and bending wire up to .030"/0.76mm.



**DE-1963**  
Loop forming in four length 3, 4, 5 & 6 mm wire capacity 0.22", 56mm.



**DE-1957**  
Johnson pliers for crown band contouring



**DE-1964**  
Arch bending plier design for convenient handling of square or rectangular wire, .020" / 0.51mm.






# PLIERS FOR ORTHODONTICS & PROSTHETICS



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروتستیس


## LAP & L-KEY JOINT PLIERS



**DE-1965**  
In order to skirt  
bands and crowns.



**DE-1966**  
Wires up to  
.020" / 0.7mm.



**DE-1967**  
Round and concave plier  
for loop forming and banding  
wire up to .20" / .051mm.



**DE-1968**  
Bending for round square  
or rectangular wire up  
to 0.30"/0.76mm



**DE-1969**  
for double back and triple  
back bands in 0.30" / .76mm  
& 0.36" / .91mm lingual arch wires.



**DE-1970**  
Detailing step plier double-sided  
offset beak permits bayonet bends.  
Creates right hand or left hand  
step-up or stepdown at 1/2, 3/4, or  
1mm increments.



**DE-1971**  
"V" bend stop special design makes  
1mm "v" bends to shorten archwires  
or provide a positive stop. Excellent  
for placing stops in titanium wires.  
For round or rectangular wires up to .025"



**DE-1972**  
Angle how plier efficient utility plier  
with serrated beaks for a firm grip.  
45 angled tips to access hard to  
reach areas.



**DE-1973**  
Weingardt plier serrated tips hold  
wires firmly at a convenient working  
angle. Tapered beaks fit easily  
between brackets and hard to reach  
areas. Our most popular utility plier.



**DE-1974**  
Pin & Ligature Cutters  
TC.



**DE-1975**  
Mini Pin & Ligature  
Cutter TC.



**DE-1976**  
Standard Ligature  
Cutter TC.



**DE-1977**  
V Bend Stop Plier  
TC.



**DE-1978**  
Angled Direct Bond Bracket  
Remover TC.




**DE-1979**  
Cutter  
Str.



Distal Cutter  
14cm  
**DE-1980** 2.5 mm  
**DE-1981** 3.0 mm  
**DE-1982** 3.0 mm



**DE-1983**  
Arrow Clasp  
Forming Pliers  
12.0cm



**DE-1984**  
Schwarz  
12.5cm

**DE-1985**  
Distal end Cutter, Safety Hold Plier  
(Ideal for Cutting Wire) TC.



**DE-1986**  
Hook Crimping Plier Angled  
Straight TC.



**DE-1987**  
3 Prong Clasp Adjustable  
Plier TC.



**DE-1988**  
Utility Plier Weingart TC.



**DE-1989**  
Direct Bracket Remover  
TC.



**DE-1990**  
Straight Direct Bond Bracket  
Remover TC.



**DE-1991**  
Pin Holding Plier TC.



**DE-1992**  
Hard wire cutters



**DE-1993**  
Elastic Separating Plier  
Round Handle



**DE-1994**  
Elastic  
Separating Plier



**DE-1995 Str.**  
**DE-1996 Cvd.**  
Smaha  
Elastomeric  
ligature  
placing forceps  
140 mm



**DE-1997**  
Elastomeric  
ligature placing  
forceps with  
groove  
Falcon 140 mm



**DE-1998**  
Elastomeric  
ligature placing  
forceps with  
groove  
Falcon 140 mm



- Щипцы для ортодонтии и протезирования
- Pinces pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروتستيس

## CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS



**DE-1999**  
Weingart Pliers



**DE-2004**  
Adhesive  
Removing Pliers



**DE-2000**  
Slim Weingart  
Pliers



**DE-2005**  
Posterior Band Removing  
Pliers - Long



**DE-2001**  
How Pliers  
Straight



**DE-2006**  
Posterior Band Removing  
Pliers - Short



**DE-2002**  
How Pliers Offset



**DE-2007**  
Bird Beak Pliers



**DE-2003**  
Band seating Pliers



**DE-2008**  
Omega Loop  
Forming Pliers

Replacement tip and  
hex key included.



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومبشیتس

## CUTTERS



**DE-2009**  
Hard  
Wire Cutter 15



**DE-2013**  
Mini Pin &  
Ligature Cutter



**DE-2010**  
Distal End Cutter  
Flush Cut



**DE-2014**  
Hard Wire Cutter  
Straight



**DE-2011**  
Universal Cut & Hold  
Distal End Cutter



**DE-2015**  
Micro-Mini Pin &  
Ligature Cutter



**DE-2012**  
Pin & Ligature  
Cutter



**DE-2016**  
Cutter

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Щипці для ортодонції та протезування
- Pinces pour l'orthodontie et la prothèse
- Alicates para ortodoncia y prótesis
- Zangen für Kieferorthopädie und Prothetik
- كماشة لتقويم الأسنان بروتستيس



**DE-2017**  
Three Jaw Pliers



**DE-2021**  
Jarabak Pliers



**DE-2018**  
Light Wore Bird Beak  
Pliers with Groove



**DE-2022**  
Hollow Chop  
Contouring Pliers



**DE-2019**  
Tweed loop  
Forming Pliers



**DE-2023**  
Lingual Arch  
Forming Pliers




**DE-2020**  
Arch Bending Pliers



**DE-2024**  
Rectangular Arch  
Bending Pliers

## ORTHODONTICS PLIERS STANDS

- Ортодонтия Щипцы Подставки
- Supports de pinces d'orthodontie
- Kieferorthopädie Zangenständer
- Ортодонтия Щипці Підставки
- KFO Zangen Ständer
- تقويم الأسنان كماشة تقف



**DE-2025**  
Plier Stand  
Stainless Steel.



**DE-2026**  
Plier Stand  
Stainless Steel.



**DE-2027**  
Pliers Stand  
Stainless Rod  
(Aluminum Anodized)  
Stainless Steel



**DE-2028**  
Plier Stand  
Stainless Steel.



**DE-2029**  
Plier Stand  
Stainless Steel.





**DE-2030**  
Plier Stand  
Plastic



## Normal Anatomy – Dental Implants

Bone Level Normal  
(No Bone Loss)



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

## IMPLANT GAGES



**DE-2100**  
**Multifunctional Meter**  
Length 90 mm,  
different sizes,  
8-10-11mm,  
corrosion-resistant

**DE-2101**  
**Multifunctional Meter**  
Length 90 mm,  
different sizes,  
5-13-15-17 mm,  
corrosion-resistant

**DE-2102**  
**Measuring Probe**  
Length 170 mm,  
hand shape or flexible  
(6-8-10-14-16 mm),  
corrosion-resistant

**DE-2103**

**DE-2104**

**DE-2105**

**DE-2106**

## IMPLANT GAGES



**DE-2107**

**DE-2108**

**DE-2109**

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantologia
- زراعة الأسنان



**DE-2110**  
Ridge Mapping  
Double Caliper



**DE-2111**  
Bone Thickness Gauge  
40mm, 180mm



**DE-2112**  
Osteometer  
Measuring Forceps For Jaw  
Ridge Width measurement, scale  
on both sides, rustfree, autoclavable




**DE-2113**  
Micro Boley  
Gauge Twisted



**DE-2114**  
Micro Boley  
Gauge Straight



**DE-2115**  
Ortho Ruler / Cone Socket




**DE-2116**  
Dynamometer




- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

## BLACK TITANIUM COLOUR COATED INSTRUMENTS



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

## BLACK TITANIUM COLOUR COATED INSTRUMENTS



**DE-2125**  
A  
Cosmetic  
Contouring

**DE-2126**  
B  
Cosmetic  
Contouring

**DE-2127**  
C  
Cosmetic  
Contouring

**DE-2128**  
D  
Cosmetic  
Contouring

**DE-2129**  
F  
Cosmetic  
Contouring

**DE-2130**  
Marginal Ridge  
& Embrasure  
Shaping Instrument

**DE-2131**  
Occlusal  
Anatomy  
Instrument

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantologia
- زراعة الأسنان

## BLACK TITANIUM COLOUR COATED INSTRUMENTS

### BURNISHERS

### FREEDMAN BURNISHERS



DE-2132  
21B



DE-2133  
27/29



DE-2134  
2  
Ladmore



DE-2135  
3  
Ladmore



DE-2136  
Freedman  
"Duckhead"



DE-2137  
Freedman  
Small Contact  
Forming



DE-2138  
Freedman  
Large Contact  
Forming




DE-2139  
Small/  
Medium Ball  
Burnisher



- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان

## SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM



DE-2140  
Sinuslift 1

DE-2141  
Sinuslift 2

DE-2142  
Sinuslift 3

DE-2143  
Sinuslift 4

DE-2144  
Sinuslift 5

DE-2145  
Sinuslift 6

DE-2146  
Membrane  
3Positioning 1

DE-2147  
Condenser 1


- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

**Sinus Lift Instrument,**  
Flexible Sinus Lift instruments  
reduce the risk of perforation of the  
"Schneider's membrane"

**Sinus Lift Elevator,**  
10 mm Straight grip  
Length 175 mm, Sinus lift elevators, dull,  
for gentle peeling and lifting the  
Schneider membrane

**Sinus lift instrument**  
10 mm Straight grip Length 180 mm,  
Buser (B1, B2, B3) raspatories,  
quality medical stainless steel,  
corrosion-resistant



DE-2148  
Fig. 1

DE-2149  
Fig. 2

DE-2150  
Fig. 3

DE-2151  
Fig. 1

DE-2152  
Fig. 2

DE-2153  
Fig. 3

DE-2154  
B01

DE-2155  
B02

DE-2156  
B03

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان

**Sinus lift instrument** 10 mm Straight grip Sinus lift elevators, blunt or sharply curved (H01, H02, K03, K04, K78, K77)



**DE-2157**  
Sinus lift  
instrument  
H01

**DE-2158**  
Sinus lift  
instrument  
H02

**DE-2159**  
Sinus lift  
instrument  
K03

**DE-2160**  
Sinus lift  
instrument  
K04

**DE-2161**  
Sinus lift  
instrument  
K78

**DE-2162**  
Sinus lift  
instrument  
K77



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

## Composite Filling Instruments Set of 6 Pcs



DE-2163



DE-2164



DE-2165



DE-2166



DE-2167



DE-2168


### COMPOSITE INSTRUMENTS USE:

A composite (tooth colored) filling is used to repair a tooth that is affected by decay, cracks, fractures, etc. The decayed or affected portion of the tooth will be removed and then filled with a composite filling.

There are many types of filling materials available, each with their own advantages and disadvantages. You and your dentist can discuss the best options for restoring your teeth. Composite fillings, along with silver amalgam fillings, are the most widely used today. Because composite fillings are tooth colored, they can be closely matched to the color of existing teeth, and are more aesthetically suited for use in front teeth or the more visible areas of the teeth.

As with most dental restorations, composite fillings are not permanent and may someday have to be replaced.

They are very durable, and will last many years, giving you a long lasting, beautiful smile.



DE-2169

## IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

Set of 5 Composite Gold Tip Instruments  
**DE-2170**



Set of 4 Composite Blue Tip Instruments  
**DE-2171**



Set of 3 Composite Blue Tip Instruments  
**DE-2172**



Set of 2 Composite Blue Tip Instruments  
**DE-2173**



**DE-2174** Periosteal Elevators Molt # 9

**DE-2175** Endodontic explorer DG # 16

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantologia
- زراعة الأسنان

## VISTA TUNNELING INSTRUMENTS SET



DE-2176



DE-2177



DE-2178



DE-2179



DE-2180



DE-2181



DE-2182





## IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie


- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



**DE-2183**  
used for elevation of sinus membrane



**DE-2184**  
used for elevation of sinus membrane



**DE-2185**  
Cobra in / out (2-4-5-6-8-10mm)



**DE-2186**  
Cobra 3.5mm



**DE-2187**  
sinus 4.8 / 5.8 (2-4-5-6-8-10mm)



**DE-2188**  
sinus 2.8 / 3.8 (2-4-5-6-8-10mm)



**DE-2189**  
Bone Graft 3.8SP / 2.8mm (2-4-5-6-8-10mm)



**DE-2190**  
Cobra in / out (2-4-5-6-8-10mm)



**DE-2191**  
Cobra 3.5mm



**DE-2192**  
Gold Color Set Dental Oral Implantology Sinus Set of 10 Pcs with cassette

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان



**DE-2193**  
Sinus Lift elevator 1



**DE-2194**  
Sinus Lift elevator 2



**DE-2195**  
Sinus Lift elevator 3



**DE-2196**  
Sinus Lift elevator 4



**DE-2197**  
Sinus Lift elevator 5



**DE-2198**  
Sinus Lift elevator 6



**DE-2199**  
Sinus Lift Elevators Set of 7 Pcs  
with cassette

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantologia
- زراعة الأسنان


## Sinus Left Elevator Set of 10 Pcs





- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

## Blue Titanium Coated Bone Grafting Instruments



**DE-2216**  
Blue Titanium Coated Bone Grafting Instruments Set of 7 Pcs with cassette

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantologia
- زراعة الأسنان



DE-2217  
Fig. 1



DE-2218  
Fig. 2



DE-2219  
Fig. 3



DE-2220  
Fig. 4



DE-2221  
Fig. 5



DE-2222  
Fig. 6



DE-2223



DE-2224




DE-2225



DE-2226



DE-2227



## Flexible Sinus-Lift Curettes – for an optimum sensitivity during the direct lift procedure

The perforation of the Schneiderian membrane is a complication that often occurs during a maxillary sinus floor augmentation procedure. The Sinus-Lift Instruments increase the sensitivity by blunt and flexible working tips with "Swing-Effect". Because of their special shapes and the characteristic of the stainless steel between the handle and the active part, light pressure is cushioned and gives the operator a direct and sensitive response. At the same time the risk to perforate the Schneiderian membrane is reduced to a minimum. Already at the development the focus was a practicable handling. Because of varying angulations and shapes of the active parts also in difficult areas they are easy to navigate also in difficult areas. The line of curettes consists of six different figures for every situation. They are complemented by spoon and condenser to place the bone graft.



The "Schneider's Membrane" of the sinus can be detached with curette by 1–2mm along the total cortical bone.



Curette is particularly useful in the lower part for detaching the membrane further. The same goes for the caudal area, especially when the membrane is 4–5 mm above the base of the maxillary sinus and the maxillary sinus is built up with granular material.



Curette makes it possible to lift the membrane further in mesial region. One working tip is especially designed for the right maxillary sinus and the other for the left maxillary sinus.



With curette the sinus membrane is further elevated in the upper and posterior part of the maxillary sinus wall.



With curette the sinus membrane is further detached over a broad area in the inner region, especially along the roof of the palate.

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان

## IMPLANT SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM



DE-2228  
Acute 6577

DE-2229  
Sinus1

DE-2230  
Sinus2

DE-2231  
Sinus  
Obtuse 6577

DE-2232  
Sinus 4

DE-2233  
Sinus 5

DE-2234  
Sinus 6

DE-2235  
Sinus 7

## SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM



DE-2236  
Sinuslift 2

DE-2237  
Sinuslift 3

DE-2238  
Sinuslift 4

DE-2239  
Sinuslift

DE-2240  
Sinuslift

DE-2243  
Membrane  
3Positioning 1

DE-2244  
Condenser 1

DE-2237  
Sinuslift 1



9 Pcs of Set  
DE-2236



9 Pcs of Set  
DE-2245



# IMPLANTOLOGY



- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

IMPLANT SINUS LIFT INSTRUMENTS  
PERIODONTAL SURGERY SURGICAL CHISELS  
KRAMER NEVINS SINUS LIFT

## SURGICAL CHISELS



DE-2246  
TG



DE-2247 SOLT3/4



DE-2248 TGO



DE-2249 OCH3



DE-2250 OCH4



DE-2251 FEDI1



DE-2252 FEDI2



DE-2253 FEDI3



DE-2254 13K/TG



DE-2255  
13K/TG Titanium



DE-2256 13K/TG



DE-2257  
13KL/R Titanium




DE-2258 13KL 4/6mm  
Titanium Available



DE-2259  
Cassette For Instruments

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantologia
- زراعة الأسنان



FLEXIBLE PERITOME

Soft Brushing Kit Set of 6 Pcs



# IMPLANTOLOGY



- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان

## SINUS ELEVATOR & MEMBRANE LIFT TIPS



**DE-2265**  
Membrane Detach Elevator



**DE-2266**  
Membrane Elevator



**DE-2267**  
Membrane Elevator



**DE-2268**  
Membrane Elevator



**DE-2269**  
Membrane Elevator



**DE-2270**  
Bone Graft Packer



**DE-2270A** Pt2 Titanium  
A/V Regular



**DE-2270B** Pt1 Titanium  
A/V Regular



**DE-2270C** Pt3 Titanium Angel S/E



**DE-2270D** Pt4 Titanium STR S/E



**DE-2270E** Pt5 Titanium TG CVD S/E



**DE-2271**  
Sinus Left Elevator & Membrane Lift Tips 3  
set of 7 Pcs



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

### MIXING SPATULA WITH SHOVEL



**DE-2272**  
Bone Carrier  
Titanium



**DE-2273**  
Raspator, Prichard  
Length 180 mm,  
10 mm hollow handle,  
double-ended  
10.5 mm | 4.5 mm,



**DE-2274**  
Raspator, Prichard  
Length 180 mm,  
10 mm hollow handle,  
double ended



**DE-2275**  
Raspator, Prichard  
Length 180 mm,  
10 mm hollow handle,  
double-ended  
10.5 mm | 4.5 mm,



**DE-2276**  
Mixing Spatula  
with shovel  
17cm



**DE-2277**  
Mixing Spatula  
with shovel  
18.5cm



**DE-2278**  
Graft Material Spoon  
Concave Tip For  
Larger Grafts

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

## BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS

Product offering includes instruments for harvesting, mixing and compacting bone Provides clinicians with the proper instrumentation for bone grafting procedures Weighted bone basin Prevents tipping and spilling of bone and other materials during mixing



**DE-2279**

Graft Packer, Concave Tips DE 4/6mm



**DE-2280**

Bone Spoon/ Packer DE 3mm



**DE-2281**

Blade and pusher Length 180 mm, for straight or angled, blade: 15 x 5 mm, pusher: 4 mm diameter Straight



**DE-2282**

Blade and pusher Length 180 mm, for straight or angled, blade: 15 x 5 mm, pusher: 4 mm diameter angular



**DE-2283**

Spoon and pusher Length 180 mm, for angled spoon: 10 x 5 mm, pusher: diameter 3 mm



**DE-2284**

Graft Packer, Smooth Tips DE 5/8mm



**DE-2285**

Graft Packer, Serrated Tips DE 3/4mm



**DE-2286**

Graft Packer, Serrated Tips DE 6/8mm



**DE-2287**

Titanium graft packer/placement site depth gauge 3.3mm & 4.0mm Diameter tips Marked at 8, 11, 13, 15 & 18mm




**DE-2288**

Tamper Length 170 mm, for bone material, double-ended, diameter 3.2 mm (2-4-6-8-10) and diameter 2.5mm (2-4-6-8-10)

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

## SPOON AND PUSHER (TRINOVO)

Trinovo handle, different workers, for bone material, high-grade medical stainless steel, corrosion-resistant





# IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان

**Hollow Handle 10mm Length 170mm**



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

## TISSUE MEMBRANE FORCEPS



**DE-2305**  
Tissue Graft Forcep  
Serrated Teeth



**DE-2306**  
Self Locking  
Grooved



**DE-2307**  
Tissue Graft Forcep  
Serrated Teeth



**DE-2308**  
Titanium  
Locking Plier



**DE-2309**  
TC Gripper Carbide  
Tipped Tissue Forceps  
Serrated Carbide Tips  
14.5cm



**DE-2310**  
Adson Carbide  
Tipped Tissue Forceps  
Serrated Carbide Tips  
12cm



**DE-2311**  
Micor Adson 1x2  
Tissue Forcep 12.1cm  
1x2 Teeth - 0.5mm Tips



**DE-2312**  
Adson 7x7  
Tissue Forceps  
7x7 Teeth 12cm



**DE-2313**  
Adson 7x7 Curved  
Tissue Forceps  
7x7 Teeth 12cm



**DE-2314**  
Suture Pliers  
Suture Forceps with  
1.6mm Ring Tip



**DE-2315**  
Debaquey Atraumatic  
Tissue Forceps Straight  
Serrated Teeth  
16cm



**DE-2316**  
Debaquey Atraumatic  
Tissue Forceps Curved  
Serrated Teeth  
16cm

- Инструменты для имплантологии
- Implantologie
- Implantsologie

- Інструменти для імплантології
- Implantologia
- زراعة الأسنان

## TISSUE MEMBRANE FORCEPS



**DE-2317**  
Russian Long  
Tissue Forceps  
Serrated Tip  
20cm



**DE-2318**  
Russian  
Tissue Forceps  
Serrated Tip  
15cm



**DE-2319**  
Tissue Punch Forceps  
To Remove Tissue  
After Using Punch



**DE-2320**  
13.5 cm  
narrow



**DE-2321**  
10.5 cm  
**Membrane tweezers KUEHNE**



**DE-2322**  
13.5 cm  
wide



**DE-2323**  
**Pin-Tweezers**  
with or without lock,  
length 15 cm,  
with lock



**DE-2324**  
without lock



**DE-2325**  
**Haltepinzette with lock**  
15,5 cm, for pins, towards  
the tips, etc



**DE-2326**  
Carry Forceps  
15cm



**DE-2327**  
Carry Forceps  
15cm



**DE-2328**  
Carry Forceps  
15cm



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantoloia
- زراعة الأسنان


## OSTEOTOME CONVEXES / CONCAVE / CONIQUES

CONVEXES



**DE-2329**  
SET OF 5 PIECES

CONCAVES







**DE-2330**  
SET OF 5 PIECES

CONIQUES



**DE-2331**  
SET OF 5 PIECES



	CONVEXES STRAIGHT	DE-2329	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES STRAIGHT	DE-2330	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES STRAIGHT	DE-2331	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm

# IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantoloia
- زراعة الأسنان

## OSTEOTOME CONVEXES / CONCAVE / CONIQUES

CONVEXES



**DE-2332**  
**SET OF 5 PIECES**

CONCAVES





**DE-2333**  
**SET OF 5 PIECES**

CONIQUES



**DE-2334**  
**SET OF 5 PIECES**



CONVEXES CURVED	DE-2332	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
CONCAVES CURVED	DE-2333	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
CONIQUES CURVED	DE-2334	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantologia
- زراعة الأسنان

## BONE SCRAPER STRAIGHT & CURVED



**DE-2335**  
Bone Scraper  
Straight



**DE-2336**  
Bone Scraper  
Curved

## TISSUE PUNCH STRAIGHT & CURVED



**Tissue Punch**  
Straight

DE-2337  
2.0mm

DE-2339  
5.0mm

DE-2341  
7.0mm

DE-2338  
4.0mm

DE-2340  
6.0mm

DE-2342  
8.0mm



**Tissue Punch**  
Curved

DE-2343  
2.0mm

DE-2345  
5.0mm

DE-2347  
7.0mm

DE-2344  
4.0mm

DE-2346  
6.0mm

DE-2348  
8.0mm

Ø4.0 Ø5.0 Ø6.0 Ø7.0



**DE-2349**  
Tissue Punch with pin  
Set of Four Pieces



**DE-2350**  
Tissue Punch  
Set of Five Pieces



**DE-2351**  
Tissue Punch  
Set of Five Pieces



- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

## BONE MILL BONE WELL



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для імплантології
- Implantoloia
- زراعة الأسنان

## TREPHINE

Burs with depth markings are suited for precise and simple transplantation of cylindrical bone grafts as well as for all other bone-extraction tasks.



**DE-2367**  
Trephine Drill  
Set of Eight Pieces



**DE-2368**  
Trephine Drill  
Set of Eight Pieces

## BONE INJECTORS, COLLECTORS & SUCTION TUBE



# IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для імплантології
- Implantologia
- زراعة الأسنان

## BONE EXPANDER SCREW & SURGICAL KITS



Ø2.0mm Ø2.5mm Ø2.8mm

**DE-2377** Short Implant Drills



Ø3.2mm Ø3.65mm Ø4.2mm Ø5.2mm

**DE-2378** Long Implant Drills



**DE-2379**

### Bone Universal Kit

Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm  
Trephines Bit: 3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.  
Driver: 2.5mm, 2.5mm Hand Pcs, 1.25mm, 1.25mm Hand Pcs  
Short Driver 2.5mm, Long Driver 2.5mm  
1.25mm, Short Hex Driver Hand Pcs  
1.25mm, Long Hex Driver Hand Pcs, Torque Wrench 4.0mm



**DE-2380**

### Bone Expander Kit

Expander: 2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm  
Ratchet Driver, Parallel Pin  
Pilot Drill, Drill 1.4mm, Drill Driver 4mm  
Hand Driver 4mm, Ratchet 4.0mm



**DE-2381**

### Bone Universal Kit

Drill: 2mm, 2.8mm, 3.2mm, 3.65mm, 4.2mm, 4.5mm, 4.8mm, 5.2mm  
Trephines Bit: 3.0mm, 4.0mm, 5.0mm, 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm.

### Bone Compression:

2mm, 3.3mm, 2.8mm, 4.2mm, 3mm, 3.5mm  
Driver: 2.5mm, 2.5mm Hand Pcs  
1.25mm, 1.25mm Hand Pcs  
Short Driver 2.5mm  
Long Driver 2.5mm  
1.25mm, Short Hex Driver Hand Pcs  
1.25mm, Long Hex Driver Hand Pcs  
Torque Wrench 6.35mm 1 Pcs

**DE-2382**

### Conical Drills: 20 Pcs

Tissue Punches: 3.0mm, 3.50mm, 4.0mm, 4.5mm, 5.0mm, 6.0mm, 6.50mm  
Bone Expander: 1.6/3.8mm, 1.9/3.0mm, 2.1/3.6mm, 2.6/3.8mm, 2.9/4.3mm  
Drill Adapter: 1 Pcs  
1.25mm, 2.50mm Hex Driver Hand Pcs



### Pilot Drill: 1.6mm

Drill Bit: 1.4mm, 1.6mm, 1.8mm, 2.0mm, 2.2mm  
Ratchet Driver: 4mm/6.35, 1.25/2.50  
Rimer: 3.2/4.7mm, 3.7/4.7mm  
Parallel Pin: 4.5mm, 2.5mm  
Double step Drill: 2.4/2.8mm  
Depth Gauge: 1 Pcs  
Torque Wrench 6.35mm Universal: 1 Pcs



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

## KEYS FOR IMPLANT SCREW



DE-2383



DE-2384



DE-2385



DE-2386

## TORQUE WRENCH & RATCHET



DE-2387 Star Torque Wrench 6.35mm Universal



DE-2388 Round Torque Wrench 4.0mm Universal



DE-2389 Hex Driver With 4.0mm Driver



DE-2390 Star Torque Wrench 6.35mm Universal



DE-2391 Hex Torque Wrench 4.00mm Universal



DE-2392 Star Torque Wrench 6.35mm Universal



DE-2393 Star Torque Wrench 6.35mm Universal



DE-2394 Star Torque Wrench 6.35mm Universal



DE-2395 Star Torque Wrench 6.35mm Universal



DE-2396 Star Torque Wrench 6.35mm Universal



DE-2397 Star Ratchet Wrench 6.35mm Universal



DE-2398 Star Ratchet Wrench 6.35mm Universal



DE-2399 Star Ratchet Wrench 6.35mm Universal



DE-2400 Star Ratchet Wrench 4.0mm Universal



DE-2401 Round Ratchet Wrench 4.0mm Universal



DE-2402 Round Ratchet Wrench 3.35mm Universal



DE-2403 Star Ratchet Wrench 4.0mm Universal



DE-2404 Star Ratchet Wrench 6.35mm Universal

# IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для імплантології
- Implantologia
- زراعة الأسنان



**DE-2405**  
Suction Cannula Kit  
with 3 Needles



Suction Pipe



Suction Needle  
1mm



Suction Needle  
1.5mm



Suction Needle  
2mm



**DE-2406**  
17.5cm

**DE-2407**  
18cm

**DE-2408**  
coarse  
Ø 5mm

Suction tubes to collect  
the extracted  
bone material  
180mm



**DE-2409**  
fine  
Ø 9mm



**DE-2410**  
Suction  
Tubes  
180mm

**DE-2411**  
Suction  
Tubes  
180mm

**DE-2412**  
Suction  
Tubes  
180mm



**DE-2413**  
Fig. 6FR  
Frazier  
15 cm - Ø 2 mm


**DE-2414**  
Ø 3 mm  
length 180mm  
**DE-2415**  
Ø 4 mm  
length 180mm  
**DE-2416**  
Ø 5 mm  
length 180mm

## IMPRESSION TRAYS

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

### PERFORATED WITH RIM LOCK

Set of 06 Pcs



#### DE-2417 SET OF 06 PCS

Upper: S, M, L  
Lower: S, M, L

### PERFORATED WITH RIM LOCK

Set of 8 Pcs



#### DE-2418 SET OF 08 PCS

Upper: XS, S, M, L  
Lower: XS, S, M, L

#### DE-2419 SET OF 08 PCS

Upper: S, M, L, XL  
Lower: S, M, L, XL

#### DE-2420 SET OF 10 PCS

Upper: XS, S, M, L, XL  
Lower: XS, S, M, L, XL



## IMPRESSION TRAYS

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel

- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

### SOLID WITH RIM LOCK

Set of 6 Pcs



#### DE-2421 SET OF 06 PCS

Upper: S, M, L  
Lower: S, M, L

### SOLID WITH RIM LOCK

Set of 8 Pcs



#### DE-2422 SET OF 08 PCS

Upper: XS, S, M, L  
Lower: XS, S, M, L

#### DE-2423 SET OF 08 PCS

Upper: S, M, L, XL  
Lower: S, M, L, XL

#### DE-2424 SET OF 10 PCS

Upper: XS, S, M, L, XL  
Lower: XS, S, M, L, XL

## IMPRESSION TRAYS

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

### PERFORATED FULL DENTURE

Set of 8 Pcs



**DE-2425 SET OF 08 PCS**  
Upper: S, M, L, XL  
Lower: S, M, L, XL

**DE-2426 SET OF 08 PCS**  
Upper: S, M, L  
Lower: S, M, L

### SOLID FULL DENTURE

Set of 8 Pcs




**DE-2427 SET OF 08 PCS**  
Upper: S, M, L, XL  
Lower: S, M, L, XL

**DE-2428 SET OF 06 PCS**  
Upper: S, M, L  
Lower: S, M, L

## IMPRESSION TRAYS

- Оттисные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع

### Perforated Without Rim Lock



### Solid With Rim Lock







I N S T R U M E N T S





**Scalers Tray**  
DE-2500 200 x 100 mm



**Scalers Trays**  
DE-2501 205 x 115 x 17 mm  
DE-2502 205 x 152 x 17 mm  
DE-2503 290 x 115 x 17 mm



**Instruments Trays**  
DE-2504 200 x 110 x 30 mm  
DE-2505 215 x 115 x 35 mm  
DE-2506 220 x 130 x 30 mm  
DE-2507 230 x 125 x 35 mm  
DE-2508 245 x 140 x 40 mm  
DE-2509 285 x 160 x 45 mm  
DE-2510 315 x 175 x 50 mm  
DE-2511 330 x 185 x 50 mm



**Instruments Trays**  
DE-2512 150 x 100 x 50 mm  
DE-2513 200 x 100 x 50 mm  
DE-2514 220 x 120 x 50 mm  
DE-2515 240 x 140 x 50 mm  
DE-2516 200 x 150 x 50 mm  
DE-2517 250 x 200 x 50 mm  
DE-2518 300 x 250 x 50 mm  
DE-2519 300 x 200 x 50 mm  
DE-2520 350 x 250 x 50 mm  
DE-2521 400 x 300 x 50 mm  
DE-2522 450 x 300 x 50 mm  
DE-2523 450 x 350 x 50 mm



**Instruments Trays With Lid**  
DE-2524 150 x 100 x 50 mm  
DE-2525 200 x 100 x 50 mm  
DE-2526 220 x 120 x 50 mm  
DE-2527 240 x 140 x 50 mm  
DE-2528 200 x 150 x 50 mm  
DE-2529 250 x 200 x 50 mm  
DE-2530 300 x 250 x 50 mm  
DE-2531 300 x 200 x 50 mm  
DE-2532 350 x 250 x 50 mm  
DE-2533 400 x 300 x 50 mm  
DE-2534 450 x 300 x 50 mm  
DE-2535 450 x 350 x 50 mm



**Mayo Tray Perforated**  
DE-2536 150 x 100 x 25 mm  
DE-2537 200 x 100 x 25 mm  
DE-2538 220 x 120 x 25 mm  
DE-2539 240 x 140 x 25 mm  
DE-2540 200 x 150 x 25 mm  
DE-2541 250 x 200 x 25 mm  
DE-2542 300 x 250 x 25 mm  
DE-2543 300 x 200 x 25 mm  
DE-2544 350 x 250 x 25 mm  
DE-2545 400 x 300 x 25 mm  
DE-2546 450 x 300 x 25 mm  
DE-2547 450 x 350 x 25 mm



**Scaler Tray**  
DE-2548 200 x 115 mm  
DE-2549 205 x 150 mm  
DE-2550 290 x 150 mm  
DE-2551 300 x 200 mm



**Instruments Tray Perforated**  
DE-2552 150 x 100 x 50 mm  
DE-2553 200 x 100 x 50 mm  
DE-2554 220 x 120 x 50 mm  
DE-2555 240 x 140 x 50 mm  
DE-2556 200 x 150 x 50 mm  
DE-2557 250 x 200 x 50 mm  
DE-2558 300 x 250 x 50 mm  
DE-2559 300 x 200 x 50 mm  
DE-2560 350 x 250 x 50 mm  
DE-2561 400 x 300 x 50 mm  
DE-2562 450 x 300 x 50 mm  
DE-2563 450 x 350 x 50 mm

## Sterilizing boxes & instruments Trays



**Instruments Box with Hooks**

DE-2564	170 x 70 x 30 mm
DE-2565	180 x 80 x 40 mm
DE-2566	200 x 100 x 40 mm
DE-2567	200 x 100 x 50 mm
DE-2568	220 x 120 x 50 mm
DE-2569	250 x 150 x 50 mm
DE-2570	250 x 200 x 50 mm
DE-2571	300 x 150 x 60 mm
DE-2572	300 x 200 x 50 mm
DE-2573	350 x 250 x 50 mm
DE-2574	400 x 200 x 100 mm
DE-2575	450 x 250 x 50 mm



**Instruments Box**

DE-2576	170 x 70 x 30 mm
DE-2577	180 x 80 x 40 mm
DE-2578	200 x 100 x 40 mm
DE-2579	200 x 100 x 50 mm
DE-2580	220 x 120 x 50 mm
DE-2581	250 x 150 x 50 mm
DE-2582	250 x 200 x 50 mm
DE-2583	300 x 150 x 60 mm
DE-2584	300 x 200 x 50 mm
DE-2585	350 x 250 x 50 mm
DE-2586	400 x 200 x 100 mm
DE-2587	450 x 250 x 50 mm



**Instruments Box Perforated**

DE-2588	170 x 70 x 30 mm
DE-2589	180 x 80 x 40 mm
DE-2590	200 x 100 x 40 mm
DE-2591	200 x 100 x 50 mm
DE-2592	220 x 120 x 50 mm
DE-2593	250 x 150 x 50 mm
DE-2594	250 x 200 x 50 mm
DE-2595	300 x 150 x 60 mm
DE-2596	300 x 200 x 50 mm
DE-2597	350 x 250 x 50 mm
DE-2598	400 x 200 x 100 mm
DE-2599	450 x 250 x 50 mm



**Gloves Box**

DE-2600 300 x 150 mm



**Catheters Tray with lid**

DE-2601	200 x 60 x 40 mm
DE-2602	200 x 90 x 40 mm
DE-2603	320 x 95 x 50 mm
DE-2604	328 x 90 x 50 mm
DE-2605	365 x 90 x 50 mm
DE-2606	440 x 115 x 50 mm



**Box With Filter**

DE-2607	100 x 50 mm
DE-2608	250 x 150 mm
DE-2609	150 x 65 mm
DE-2610	200 x 100 mm
DE-2611	300 x 200 mm



**Instruments Trays With Cover**

DE-2612	210 x 110 x 50 mm
DE-2613	240 x 160 x 50 mm
DE-2614	300 x 180 x 50 mm
DE-2615	360 x 210 x 50 mm
DE-2616	400 x 230 x 55 mm

**Instruments Trays Without Cover**

DE-2617	210 x 110 x 50 mm
DE-2618	240 x 160 x 50 mm
DE-2619	300 x 180 x 50 mm
DE-2620	360 x 210 x 50 mm
DE-2621	400 x 230 x 55 mm



**Instruments Box with lid**

DE-2622	170 x 70 x 30 mm
DE-2623	180 x 80 x 40 mm
DE-2624	200 x 100 x 40 mm
DE-2625	200 x 100 x 50 mm
DE-2626	220 x 120 x 50 mm
DE-2627	250 x 150 x 50 mm
DE-2628	250 x 200 x 50 mm
DE-2629	300 x 150 x 60 mm
DE-2630	300 x 200 x 50 mm
DE-2631	350 x 250 x 50 mm
DE-2632	400 x 200 x 100 mm
DE-2633	450 x 250 x 50 mm





**Kidney Tray**  
 DE-2634 160 mm  
 DE-2635 200 mm  
 DE-2636 250 mm  
 DE-2637 280 mm  
 DE-2638 300 mm



**Kidney Tray Shallow**  
 DE-2639 170 mm  
 DE-2640 220 mm  
 DE-2641 240 mm



**Kidney Tray With Cover**  
 DE-2642 160 mm  
 DE-2643 200 mm  
 DE-2644 250 mm  
 DE-2645 280 mm  
 DE-2646 300 mm



**Kidney Tray**  
 DE-2647 160 mm  
 DE-2648 200 mm  
 DE-2649 250 mm  
 DE-2650 280 mm  
 DE-2651 300 mm



**Forceps Jars**  
 DE-2652 30 x 90 mm  
 DE-2653 50 x 100 mm  
 DE-2654 50 x 130 mm  
 DE-2655 50 x 175 mm  
 DE-2656 75 x 100 mm  
 DE-2657 75 x 130 mm  
 DE-2658 75 x 175 mm  
 DE-2659 100 x 175 mm  
 DE-2660 100 x 200 mm



**Methot Jar Without Lid**  
 DE-2661



**Dressing Jar / Cotton Container**  
 DE-2662 80 x 80 mm  
 DE-2663 100 x 60 mm  
 DE-2664 100 x 80 mm  
 DE-2665 100 x 100 mm  
 DE-2666 120 x 80 mm  
 DE-2667 120 x 100 mm  
 DE-2668 150 x 100 mm  
 DE-2669 150 x 150 mm  
 DE-2670 180 x 150 mm  
 DE-2671 180 x 180 mm  
 DE-2672 200 x 200 mm  
 DE-2673 240 x 240 mm



**Sputum Cup**  
 DE-2674



**Dressing Jar With Lid & Knob**  
 DE-2675 80 x 80 mm  
 DE-2676 90 x 90 mm  
 DE-2677 100 x 100 mm  
 DE-2678 120 x 120 mm  
 DE-2679 100 x 130 mm  
 DE-2680 150 x 150 mm  
 DE-2681 125 x 170 mm



**Round Bowls**  
 DE-2682 80 x 35 mm  
 DE-2683 125 x 60 mm



**Spirit Lamp**  
DE-2684



**Spirit Lamp**  
DE-2685



**Spirit Lamp**  
DE-2686



**Cotton Holding Jar**  
DE-2687



**Cotton Holder**  
DE-2688



**Cotton Holder**  
DE-2689



**Cotton Wool Dispenser Acc.  
to Methot**  
DE-2690



**Needle Sterillization Tub**  
DE-2691



**Round Endodontic Boxes**  
complete with lid  
DE-2692 Ø mm. 94 - h. mm. 70



**Endodontic Box**  
DE-2693



**DE-2694**  
100 x 180 x 11mm  
Mini size



**DE-2695**  
100 x 180 x 11mm  
Mini size



**DE-2696**  
100 x 180 x 11mm  
Mini size



**DE-2697**  
100 x 180 x 11mm  
Mini size



**DE-2698**  
100 x 180 x 11mm  
Mini size



Lid  
Coperchio

**DE-2699**  
103 x 182 x 10mm  
Mini size



**DE-2700**  
Silver



**DE-2701** Blue  
**DE-2702** Red



**DE-2703** Blue  
**DE-2704** Red



**DE-2705**  
Silver



**DE-2706**  
Green



**DE-2708**  
Blue



**DE-2707**  
Golden



**DE-2709**  
Red



**DE-2710**  
Silver



**DE-2711**  
Midi  
125 x 105 x 54mm



**DE-2712**  
Maxi  
204 x 105 x 54mm



**DE-2713**  
Maxi  
204 x 105 x 54mm





## Round Sterilizer Drum

DE-2714	mm. Ø 80 x 80 H
DE-2715	mm. Ø 125 x 80 H
DE-2716	mm. Ø 125 x 100 H
DE-2717	mm. Ø 125 x 125 H
DE-2718	mm. Ø 150 x 80 H
DE-2719	mm. Ø 150 x 100 H
DE-2720	mm. Ø 150 x 130 H
DE-2721	mm. Ø 150 x 150 H
DE-2722	mm. Ø 150 x 190 H
DE-2723	mm. Ø 190 x 125 H
DE-2724	mm. Ø 190 x 145 H
DE-2725	mm. Ø 190 x 165 H
DE-2726	mm. Ø 190 x 190 H
DE-2727	mm. Ø 190 x 240 H
DE-2728	mm. Ø 190 x 290 H
DE-2729	mm. Ø 240 x 120 H
DE-2730	mm. Ø 240 x 145 H
DE-2731	mm. Ø 240 x 165 H
DE-2732	mm. Ø 240 x 190 H
DE-2733	mm. Ø 240 x 240 H
DE-2734	mm. Ø 240 x 290 H
DE-2735	mm. Ø 240 x 340 H
DE-2736	mm. Ø 240 x 390 H
DE-2737	mm. Ø 290 x 145 H
DE-2738	mm. Ø 290 x 165 H
DE-2739	mm. Ø 290 x 190 H
DE-2740	mm. Ø 290 x 240 H
DE-2741	mm. Ø 290 x 290 H
DE-2742	mm. Ø 290 x 340 H
DE-2743	mm. Ø 290 x 390 H
DE-2744	mm. Ø 340 x 145 H
DE-2745	mm. Ø 340 x 165 H
DE-2746	mm. Ø 340 x 190 H
DE-2747	mm. Ø 340 x 240 H
DE-2748	mm. Ø 340 x 290 H
DE-2749	mm. Ø 340 x 340 H
DE-2750	mm. Ø 340 x 390 H
DE-2751	mm. Ø 390 x 140 H
DE-2752	mm. Ø 390 x 165 H
DE-2753	mm. Ø 390 x 190 H
DE-2754	mm. Ø 390 x 240 H
DE-2755	mm. Ø 390 x 290 H
DE-2756	mm. Ø 390 x 340 H
DE-2757	mm. Ø 390 x 390 H
DE-2758	mm. Ø 390 x 490 H



## Round Sterilizer Drums "Filter"

DE-2759	mm. Ø 80 x 80 H
DE-2760	mm. Ø 125 x 80 H
DE-2761	mm. Ø 125 x 100 H
DE-2762	mm. Ø 125 x 125 H
DE-2763	mm. Ø 150 x 100 H
DE-2764	mm. Ø 150 x 130 H
DE-2765	mm. Ø 150 x 150 H
DE-2766	mm. Ø 150 x 190 H
DE-2767	mm. Ø 190 x 145 H
DE-2768	mm. Ø 190 x 165 H
DE-2769	mm. Ø 190 x 190 H
DE-2770	mm. Ø 190 x 240 H
DE-2771	mm. Ø 240 x 145 H
DE-2772	mm. Ø 240 x 165 H
DE-2773	mm. Ø 240 x 190 H
DE-2774	mm. Ø 240 x 240 H
DE-2775	mm. Ø 240 x 290 H
DE-2776	mm. Ø 240 x 340 H
DE-2777	mm. Ø 240 x 390 H
DE-2778	mm. Ø 290 x 145 H
DE-2779	mm. Ø 290 x 165 H
DE-2780	mm. Ø 290 x 190 H
DE-2781	mm. Ø 290 x 240 H
DE-2782	mm. Ø 290 x 290 H
DE-2783	mm. Ø 290 x 340 H
DE-2784	mm. Ø 290 x 390 H
DE-2785	mm. Ø 340 x 145 H
DE-2786	mm. Ø 340 x 165 H
DE-2787	mm. Ø 340 x 190 H
DE-2788	mm. Ø 340 x 240 H
DE-2789	mm. Ø 340 x 290 H
DE-2790	mm. Ø 390 x 140 H
DE-2791	mm. Ø 390 x 165 H
DE-2792	mm. Ø 390 x 190 H
DE-2793	mm. Ø 390 x 240 H
DE-2794	mm. Ø 390 x 290 H
DE-2795	mm. Ø 390 x 340 H
DE-2796	mm. Ø 390 x 490 H




## Round Bowls With Cover

DE-2797	Ø 40 mm ml.30
DE-2798	Ø 60 mm ml.50
DE-2799	Ø 80 mm ml.150
DE-2800	Ø 100 mm ml.250
DE-2801	Ø 120 mm ml.450
DE-2802	Ø 150 mm ml.750
DE-2803	Ø 180 mm ml.1500
DE-2804	Ø 200 mm ml.2100
DE-2805	Ø 250 mm ml.3500

## Round Bowls With Out Cover

DE-2806	Ø 40 mm ml.30
DE-2807	Ø 60 mm ml.50
DE-2808	Ø 80 mm ml.150
DE-2809	Ø 100 mm ml.250
DE-2810	Ø 120 mm ml.450
DE-2811	Ø 150 mm ml.750
DE-2812	Ø 180 mm ml.1500
DE-2813	Ø 200 mm ml.2100
DE-2814	Ø 250 mm ml.3500

## INSTRUMENTS MANAGEMENT SYSTEM



Cassette For 4 Instruments  
DE-2815



Cassette For 5 Instruments  
DE-2816



Cassette For 6 Instruments  
DE-2817



Cassette For 7 Instruments  
DE-2818



Cassette For 10 Instruments  
DE-2819



Cassette For 20 Instruments  
DE-2820



Cassette For Extraction Forceps  
DE-2821



Cassette For Root Elevator  
DE-2822



DE-2823  
Cassette Tray  
15 Pcs



DE-2824  
Cassette Tray  
15 Pcs



## STERILIZING MESH TRAY

Different sizes & dimensions  
are available at buyer's demand



**MESH TRAYS**  
**DE-2825**

Different sizes & dimensions  
are available at buyer's demand



**MESH TRAYS**  
**DE-2826**  
**05 Pcs**



**MESH TRAYS**  
**DE-2827**  
**07 Pcs**



**MESH TRAYS**  
**DE-2828**  
**10 Pcs**



**MESH TRAYS**  
**DE-2829**  
**20 Pcs**



**DE-2830**  
Mesh Cassettes  
5 Pieces, 7 Pieces, 10 Pieces,  
15 Pieces & 20 Pieces



### MESH TRAY

DE-2831 240x130x50mm  
DE-2832 240x130x60mm  
DE-2833 240x250x50mm  
DE-2834 240x250x60mm  
DE-2835 480x250x50mm  
DE-2836 480x250x60mm  
DE-2837 480x250x70mm



### MESH TRAY

DE-2838 240x250x 50mm  
DE-2839 240x250x 60mm  
DE-2840 240x250x 70mm  
DE-2841 480x250x 50mm  
DE-2842 480x250x 60mm  
DE-2843 480x250x 70mm  
DE-2844 480x250x100mm



### MESH TRAY

DE-2845 240x250x50mm  
DE-2846 240x250x60mm  
DE-2847 240x250x70mm  
DE-2848 280x250x50mm  
DE-2849 280x250x60mm



### MESH TRAY

DE-2850 280x250x 70mm  
DE-2851 280x250x100mm  
DE-2852 280x250x 50mm  
DE-2853 280x280x 60mm  
DE-2854 280x280x 70mm  
DE-2855 300x280x 60mm





**DE-2856**  
18×7×3.5cm  
5 Instruments Cassette  
Detachable Lids  
+  
5 Instruments Cassette Hinged  
Lids with button Lock



**DE-2857**  
29×18.5×3.5cm  
20 Instruments Cassette Hinged  
Lids with surgical kit fitting,  
button lock




**Instruments Cassette**  
**DE-2858** 5 Pieces Cassette  
**DE-2859** 7 Pieces Cassette  
**DE-2860** 10 Pieces Cassette



**DE-2861**  
18×11×3.5cm  
Extra large cassette for 6  
Instruments, with compartment  
for small parts hinged lids



**DE-2862**  
29×18.5×3.5cm  
20 Instruments Cassette  
Detachable Lids  
+  
20 Instruments Cassette Hinged  
Lids with button Lock



**DE-2863**  
10×13×2.2cm  
10 Instruments Cassette Hinged  
Lids, Click Lock Available in 2  
designs



**DE-2864**  
29×18.5×3.5cm  
20 Instruments Cassette Detachable  
Lids with elevator fitting  
+  
20 Instruments Cassette Hinged  
Lids with elevator fitting , button lock



**DE-2865**  
18×9×2.2cm  
7 Instruments Cassette Hinged  
Lids, Click Lock, Available in 2  
designs



**DE-2866**  
21×14×3.5cm  
15 Instruments cassette,  
Detachable lids, button lock



**DE-2867**  
21×14×3.5cm  
15 Instruments cassette,  
detachable lids, button  
lock, with fitting of micro  
surgery kit



**DE-2868**  
21×14×4.5cm  
8 pcs luxator/elevator  
cassette, Detach0  
ablelids, button lock



**DE-2869**  
15×9×5cm  
Box for dental implant tools and  
surgical kits, Capacity:8  
tools, 18 tools, 24 Tools, 58 Tools  
+ removable compartment



**Instruments Tray Perforated**  
**DE-2870** 288 x 187 x 40 mm



**Instruments Tray Perforated Base**  
**DE-2871** 288 x 187 x 40 mm



**Instruments Tray Solid Base**  
**DE-2872** 288 x 187 x 40 mm



**DE-2873** Implant Basic Starter Kit 11 instruments with case



**DE-2874** Crown Bridge Set Instruments with Case



**DE-2875** Bone Curettes Spoon Style set of 4 with Metal Case



**DE-2876** Peritome Set of 2 with Blue Nitrate Tips and Metal Case



**DE-2877** Root Fragment Removal elevators, Set of 4 with Metal Case



**DE-2878** Peritome Set of 3 Gold Handles Serrated Blades with Metal Case



## DESCRIPTION

## PAGE#

## A

A COSMETIC CONTOURING IMPLANTOLOGY	148
AB1 BOGHOSIAN IMPLANTOLOGY	147
AB2 BOGHOSIAN IMPLANTOLOGY	147
ABOU-RASS EXPLORERS	54
ACUTE IMPLANTOLOGY	161
ADAMS PLIER	133
ADAMS PLIER ROUND BOX	133
ADAMS PLIERS FOR HOOKS AND BENDING	133
ADAMS, SMALL WIRE BENDING AND LOOP FORMING PLIERS	132
ADERER CLAMP THREE JAW WIRE BENDING PLIER	138
ADERER FACEBOW LOOP FORMING PLIERS	134
ADERER LARGE WIRE BENDING AND LOOP FORMING PLIERS	133
ADERER PLIER	134
ADHESIVE REMOVING PLIERS	36,141
ADSON 7X7 CURVED TISSUE FORCEPS	169
ADSON BROWN TWEEZERS	53
ADSON CARBIDE TIPPED TISSUE FORCEPS	169
ADSON NEEDLE HOLDERS	107
ADSON-TC TWEEZERS	53
ADULT PALATAL DENTAL PHOTOGRAPHIC MIRRORS	98
ALLIS HAEMOSTATIC FORCEPS	110
ALLIS-BABY HAEMOSTATIC FORCEPS	110
ALL-PURPOSE PLIERS	109
ALLUMINUM CALIPER	127
ALM RETRACTOR	102
AMALGAM CARRIERS	117
AMALGAM INSTRUMENTS	116,117
AMALGAM WELL	115,117
ANATOMIC HANDLE EXPLORERS	54
ANDREWS FILLING INSTRUMENTS	87
ANGLE HOW PLIER	139
ANGLED DIRECT BOND BRACKET REMOVER TC	139
ANGLED EXTRACTING FORCEPS FOR UPPER ROOTS	03
ANGLED UTILITY PLIERS	136
ANISWORTH RUBBER DAM CLAMPS FORCEPS	111
ANTERIOR	58
APEX FORCEPS LOWER ROOT DEEP ROOT PICK	06
APEX FORCEPS LOWER UNIVERSAL DEEP ROOT PICK	05
APEXO ROOT ELEVATORS	29
APICAL RETENTION FORCEPS	23
APICAL ROOT ELEVATORS	31,33
ARCH BENDING PLIERS	38,143
ARROW CLASP FORMING PLIERS	139
ARTICULATORS	87
ASPIRATING SYRING	27
ATRAUMATIC TWEEZERS	52
ATRAUMATIC WITH LOCK TWEEZERS	52

## B

BABY TOOTH EXTRACTION FORCEPS SET OF 7 PIECES	08
BACK ACTION FILLING INSTRUMENTS	88
BACKHAUS HAEMOSTATIC FORCEPS	110
BACKHAUS TOWEL CLAMPS	110
BACOLITE HANDLE WAX KNIVES	96
BALL BURNISHERS	71
BALL FILLING INSTRUMENTS	93
BALL SOCKET NOT-PERFORATING AND NOT-DAMAGING	110
BALLACORN FILLING INSTRUMENTS	88
BAND REMOVING PLIERS	32,136
BAND SEATER AUTO CLEAVABLE	129
BAND SEATER AUTO CLEAVABLE ROUND SHAPE	129
BAND SEATER AUTO CLEAVABLE SQUARE SHAPE	129
BAND SEATER AUTO CLEAVABLE TRIANGLE SHAPE	129
BAND SEATING & REMOVING INSTRUMENTS	129
BAND SEATING & REMOVING SILICON GRIP	129
BAND SEATING PLIERS	141
BANDS & CROWNS PLIER	134
BARRAQUER T.C NEEDLE HOLDERS	109
BASINS FOR BONE	174
BAUMGARTNER NEEDLE HOLDERS	06,107
BAYONET ROOTS ELEVATOR SET	36
BEAVERTAIL FILLING INSTRUMENTS	88
BEEBEE GUM SCISSORS	106
BEERENDONK CALIPER	127
BEIN ROOT ELEVATORS	29,40
BENDING PLIER	138
BESEK COMPOSITE INSTRUMENTS	72
BESEK COMPOSITE MODELLING INSTRUMENTS	85
BESEK COMPOSITE PLUGGERS INSTRUMENTS	82
BESEK PLUGGERS	76
BEYER BONE ROUNDER	25
BINANGLE SCALER	67
BIRD BEAK PLIER	38,141
BIRD BEAK PLIER FOR BENDING	138
BIRD BEAK PLIER STANDARD	133
BIRD BEAK WIRE BENDING & CUTTING PLIERS	134
BIRD BEAK WIRE BENDING PLIER	133
BISHOP RETRACTORS	101
BLACK FILLING INSTRUMENTS	80,87
BLACK SCALER	66
BLACK TITANIUM COLOUR COATED INSTRUMENTS	147,148,149
BLADE AND PUSHER	166
BLUE TITANIUM COATED BONE GRAFTING INSTRUMENTS	59,187
SET OF 7 PCS WITH CASSETTE	159
BLUMENTHAL BONE ROUNDER	25
BOLEY GAUGE	127
BONE BRACKET GAUGE	125
BONE CARRIER TITANIUM	165
BONE CRUSHER	174
BONE CRUSHER TITANIUM METAL CUTTER	174
BONE EXPANDER KIT	176
BONE EXPANDER SCREW & SURGICAL KITS	176
BONE FILES	119
BONE GRAFT	159
BONE GRAFT IMPLANTOLOGY	156
BONE GRAFT PACKER	164
BONE GRAFT/CARRIER TITANIUM	167
BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS	166
BONE INJECTORS, COLLECTORS & SUCTION TUBE	175
BONE MILL	174
BONE ROUNDER	25
BONE SCRAPER	173
BONE SCRAPER STRAIGHT & CURVED	173
BONE SPOON PACKER	166

## DESCRIPTION

## PAGE#

BONE THICKNESS GAUGE	146
BONE UNIVERSAL KIT	176
BONE WELL	174
BOX FOR DENTAL IMPLANT TOOLS AND SURGICAL KITS	191
BOX JOINT PLIERS	132,133,134
BRACKET HOLDER TWEEZERS	128
BRACKET PLACER TWEEZERS	128
BRACKET PLACING TWEEZERS	128
BRACKET REMOVING PLIER CURVED	132
BRACKET REMOVING PLIER ROUND BOX	132
BRACKET REMOVING PLIERS ANGLED	136
BRACKET REMOVING PLIERS STRAIGHT	136
BRACKET REMOVING PLIERS WITH PLASTIC TIP	136
BRACKET TWEEZERS SLOT ALIGNER	128
BREWER RUBBER DAM CLAMPS FORCEPS	111
BRUENINGS RETRACTORS	101
BUCCAL AND LINGUAL/PALATAL SCALER	61
BUCCAL DENTAL PHOTOGRAPHIC MIRRORS	98
BUCCAL SURFACES SCALER	61
BUCCAL-SMALL DENTAL PHOTOGRAPHIC MIRRORS	98
BURNISHER FILLING INSTRUMENTS	78,88
BURNISHERS	149
BURNISHERS CONSERVATIVE TREATMENT	68,71
BUSER PERIODONTAL INSTRUMENTS	121

## C

COSMETIC CONTOURING IMPLANTOLOGY	148
CARBIDE DIAMOND COATED POINT	24
CARRY FORCEPS	170
CARVER COMPOSITE MODELLING INSTRUMENTS	85
CARVERS INSTRUMENTS	89,90
CARVERS INSTRUMENTS CASSETTE FOR 5 INSTRUMENTS	90
CASSETTE FOR INSTRUMENTS	162
CASTROVIEJO CALIPERS	126
CASTROVIEJO GUM SCISSORS	105
CASTROVIEJO MICRO SCISSORS	108
CASTROVIEJO NEEDLE HOLDERS	107,109
CATHETERS TRAY WITH LID BOX WITH FILTER	184
CATRIDGE SYRING	27
CEMENT SPATULAS	95
CERI SAW	124
CHEEK RETAINER	101
CHILD PALATAL DENTAL PHOTOGRAPHIC MIRRORS	98
CHIP BLOWER	116
CHISELS, BONE FILES & SPOONE CURETTE	118
CHISELS SCALER	65,67
CHOMPET ROOT ELEVATOR	42
CINCH BACK / TIE BACK DISTAL END BENDING PLIERS	137
CLAW-LIFT ROOTS ELEVATOR	37
COARSE SUCTION TUBES	178
COBRA IMPLANTOLOGY	156
COHEN UNIVERSAL EXTRACTING FORCEPS (AMERICAN PATTERN)	17
COLLEGE FORCEPS	51
COLLEGE TWEEZER	53
COLLEGE WITH LOCK TWEEZERS	52
COLLINS TONGUE FORCEPS	110
COLUMBIA RETRACTOR	99
COMBINATION RETRACTORS	101
COMELLAS FORCEPS ANATOMICAL HANDLE	26
COMPASSO DE WILLIS	127
COMPOSITE BALL PLUGGERS	76
COMPOSITE FILLING INSTRUMENTS	92,93
COMPOSITE FILLING INSTRUMENTS SET OF 6 PCS	93
COMPOSITE GOLD TIP INSTRUMENTS	93
COMPOSITE INSTRUMENTS	72,73
COMPOSITE INSTRUMENTS BLUE TITANIUM WITH CASSETTES	89
COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC	89
COMPOSITE MODELLING INSTRUMENTS	85
COMPOSITE PLUGGERS INSTRUMENTS	82,83
COMPOSITE/PLASTIC FILLING INSTRUMENTS	87
CONDENSER IMPLANTOLOGY	150,161
CONE-SOCKET MOUTH MIRRORS	49
CONE-SOCKET FRONT SURFACE MOUTH MIRRORS	49
CONE-SOCKET MIRROR HANDLES FOR MOUTH MIRRORS	50
CONE-SOCKET PLAIN MOUTH MIRRORS	49
CONE-SOCKET RHODIUM MOUTH MIRRORS	49
CONE-SOCKET SOLID HANDLES FOR MOUTH MIRRORS	50
CONE-SOCKET SOLID OCTAGONAL HANDLES FOR MOUTH MIRRORS	50
CONTOURING PLIER	134,138
COPERCHO HOLLOWARE INSTRUMENTS	145
CORROSION-RESISTANT	186
COTTON HOLDER	186
COTTON HOLDING JAR	186
COTTON WOOL DISPENSER ACC. TO METHOT	186
COUPLAND ROOT ELEVATORS	34
CRANE KAPLAN TWEEZERS	53
CRILE-WOOD NEEDLE HOLDERS	106,107
CROWN REMOVER	114
CROWN SPREADERS	114
CRUCIBLE TONGS	124
CRYER EXTRACTING FORCEPS (AMERICAN PATTERN)	12,13,19,20
CRYER ROOT ELEVATORS	32
CUTTER	139,142
CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS	141

## D

DEAN GUM SCISSORS	105
DEBAKEY ATRAUMATIC TISSUE FORCEPS	169
DELUX Boley GAUGE	127
DELUXE BRACKET HOLDER	128
DELUXE BRACKET HOLDER TWEEZERS	128
DENHARD RETRACTORS - MOUTH GAGS	100
DENTAL PHOTOGRAPHIC MIRRORS	98
DENTAL SYRINGES	27
DENTAURUM CALIPER	127
DERF NEEDLE HOLDERS	106,107
DETAILING STEP PLIER	139
DETAILING STEP PLIERS / BAYONET BENDING PLIERS	137
DIAMOND COATED MARGIN TRIMMERS	65
DIRECT BOND BRACKET HOLDER TWEEZERS	128
DIRECT BOND BRACKET HOLDER TWEEZERS, SMALL TIP	128
DIRECT BONDING BRACKET HOLDERS	128
DIRECT BRACKET REMOVER TC	140
DISC SYSTEM DENTAL SYRINGES	27
DISSECTING SET SET OF 11 PCS	94
DISSECTING SET SET OF 9 PCS	94



DESCRIPTION	PAGE#	DESCRIPTION	PAGE#
DISTAL AMALGAM CARRIERS	117	HAMMER HEAD DISTAL END BENDING PLIERS	137
DISTAL BENDER SILICON GRIP	129	HAMMER WITH WOODEN HANDLE	124
DISTAL CUTTER	139	HANDLE ROUND HANDLES FOR MOUTH MIRRORS	50
DISTAL END CUTTER	135	HANDLES FOR MOUTH MIRRORS	50
DISTAL END CUTTER FLUSH CUT	142	HARD WIRE CUTTER	142
DISTAL END CUTTER, SAFETY HOLD PLIER (IDEAL FOR CUTTING WIRE) TC	140	HARD WIRE CUTTER STRAIGHT	142
DISTAL END, WIRE CUTTERS & BRACKET REMOVER PLIERS	140	HARD WIRE CUTTER, 15° ANGLED	135
DISTAL SCALER	61, 63	HARD WIRE CUTTER, 15° ANGLED NON-AUTOCALVABLE	135
DOUBLE ENDED CEMENT SPATULA	131	HARD WIRE CUTTER, STRAIGHT	140
DOUBLE ROUNDED JAW BENDING PLIERS	138	HARD WIRE CUTTERS	140
DOUBLE-ENDED EXPLORERS	54	HARRIS EXTRACTING FORCEPS (AMERICAN PATTERN)	13, 18
DRESSING TWEEZERS	52, 53	HARTMANN (BABY MOSQUITO) NEEDLE HOLDERS	107
DYNAMOMETER	146	HEAT CARRIERS	57
<b>E</b>		HEBER BAYONET ROOTS ELEVATOR	40
EASY ACCESS TWEEZES	128	HEBER ROOTS ELEVATOR	40
EBNETER SCALER	64	HEIDBRINK ROOT APPEX ELEVATOR	44
ELASTIC INSTRUMENT SILICON GRIP	129	HEIDBRINK ROOT ELEVATORS	33, 43
ELASTIC SEPARATING PLIER	140	HEIDBRINK ROOTS ELEVATOR WARWICK-JAMES	44
ELASTIC SEPARATING PLIER ROUND HANDLE	140	HEIDEMANN FILLING INSTRUMENTS	93
ELASTOMERIC LIGATURE PLACING FORCEPS WITH GROOVE FALCON	140	HEIDEMANN SPATULAS	73
ELEVATORS SET OF 8 PCS	46	HEIGHT GAUGE 0.18 STAINLESS STEEL	125
ELEVATORS SET OF 8 PCS	46	HEIGHT GAUGE ALUMINUM	125
ELLIOT SIGVELAND RETAINERS	115	HEMINGWAY TWEEZERS	53
ENDO FINGER RULER SPAN MEASURE SCALE	125	HEX TORQUE WRENCH UNIVERSAL	177
ENDODONTIC BOX	196	HEXAGONAL HANDLES FOR MOUTH MIRRORS	50
ENDODONTIC EXPLORER	154	HIRSCHFELD PERIODONTAL INSTRUMENTS	122
ENDODONTICS INSTRUMENTS	58	HIRT COMPOSITE MODELLING INSTRUMENTS	85
ERGO GRIP RETRACTOR	101	HIRT COMPOSITE PLUGGERS INSTRUMENTS	83
EURO BREECH LOADING DENTAL SYRINGES	27	HIRT SCALER	66
EXAVATORS EXPLORERS	51	HMG HEISTER	99
EXCAVATOR INSTRUMENTS	74	HOLLENBACK FILLING INSTRUMENTS	80, 87
EXCAVATORS	86	HOLLOW CHOP CONTOURING PLIERS	143
EXECT BRACKET HOLDER TWEEZERS	128	HOLLOW HANDLES FOR MOUTH MIRRORS	50
EXPLORERS	51, 54, 55, 56	HOLLOW MIRROR, HANDLES FOR MOUTH MIRRORS	50
EXTRACTING FORCEPS (AMERICAN PATTERN)	126, 124	HOLLOWWARE INSTRUMENTS	183, 184, 186, 187, 190, 191
EXTRACTING FORCEPS (ENGLISH PATTERN)	01, 007, 26	HOLZ SCALER	65
EXTRACTING FORCEPS (ENGLISH PATTERN) SET OF 10 PIECES	07	HOOK CRIMPING PLIER ANGLED STRAIGHT TC	140
EXTRACTING FORCEPS	08	HOOK KNIVES	96
EXTRACTING FORCEPS AMERICAN ADULT SET OF 10 PCS	28	HOPSON FILLING INSTRUMENTS	78
EXTRACTING FORCEPS AMERICAN ADULT SET OF 8 PCS	28	HOW PLIER	132
EXTRACTING FORCEPS CHILDREN PATTERN SET OF 10 PCS	28	HOW PLIER CURVED	132
EXTRACTING FORCEPS ENGLISH ADULT SET OF 10 PCS	28	HOW PLIERS OFFSET	141
EXTRACTING FORCEPS PEDO SET OF 7 PCS	28	HOW PLIERS STRAIGHT	141
EXTRACTING FORCEPS SET	28	HOW STRAIGHT UTILITY PLIERS	136
EXTRACTING FORCEPS SET OF 10 PIECES	28	HOWARD PERIOSTEAL ELEVATORS	139
EXTRACTING FORCEPS SET OF 10 PIECES	28	HULL EXTRACTING FORCEPS (AMERICAN PATTERN)	129
EXTRACTING FORCEPS SET OF 10 PIECES	28	HYLIN CEMENT SPATULAS	95
EXTRACTING FORCEPS SET OF 13 PIECES	28	HYLIN ROOTS ELEVATOR	39
EXTRA-FLEX IMPLANTOLOGY	147	<b>I</b>	
<b>F</b>		IMPLANT GAGES	145
FALCON BEAK CONTOURING PLIERS	133, 138	IMPLANT SINUS LIFT INSTRUMENTS	161, 162
FARABEUF RETRACTORS SET OF 2 PIECES	103	IMPLANTOLOGY	145, 178
FELDMAN RETRACTOR	103	IMPLANTOLOGY 9 PCS OF SET	161
FILLING INSTRUMENTS (GOLD POINTED)	91	IMPLANTOLOGY SET OF 2 COMPOSITE BLUE TIP INSTRUMENTS	154
FILLING INSTRUMENTS	87, 88, 91, 93	IMPLANTOLOGY SET OF 3 COMPOSITE BLUE TIP INSTRUMENTS	154
FILLING INSTRUMENTS CONSERVATIVE TREATMENT	69, 70, 75, 78, 79, 81, 84	IMPLANTOLOGY SET OF 4 COMPOSITE BLUE TIP INSTRUMENTS	154
FIRST AID EXAMINATION SET OF 3 PIECES	161	IMPLANTOLOGY SET OF 5 COMPOSITE GOLD TIP INSTRUMENTS	154
FLEXIBLE PERITOME	53	IMPRESSION TRAYS	179, 182
FLEX-THIN IMPLANTOLOGY	163	IMPRESSION TRAYS PERFORATED FULL DENTURE SET OF 8 PCS	181
FLOER ROOT ELEVATORS	29, 33	IMPRESSION TRAYS PERFORATED WITHOUT RIM LOCK	181
FLUSH-CUT DISTAL AND CUTTING PLIERS	133	IMPRESSION TRAYS SOLID FULL DENTURE SET OF 06 PCS	181
FLUSH-CUT DISTAL END CUTTER	135	IMPRESSION TRAYS SOLID FULL DENTURE SET OF 8 PCS	181
FOERSTER DRESSING FORCEPS	110	IMPRESSION TRAYS SOLID WITH RIM LOCK	182
FRAZIER SUCTION TUBES WITH FILTER	178	IMPRESSION TRAYS SOLID WITH RIM LOCK SET OF 06 PCS	180
FREEDMAN "DUCKHEAD"	149	IMPRESSION TRAYS SOLID WITH RIM LOCK SET OF 08 PCS	180
FREEDMAN BURNISHERS	149	IMPRESSION TRAYS SOLID WITH RIM LOCK SET OF 10 PCS	180
FREEDMAN LARGE CONTACT FORMING	149	INSTRUMENTS BOX	184
FREEDMAN SMALL CONTACT FORMING	149	INSTRUMENTS BOX PERFORATED	184
FREER PERIODONTAL INSTRUMENTS	121	INSTRUMENTS BOX WITH HOOKS	184
FREER PERIOSTEAL ELEVATORS	120	INSTRUMENTS BOX WITH LID	184
FRET SAW	124	INSTRUMENTS CASSETTE 10 HINGED	191
FRIEDMAN BONE ROUNDER	25	INSTRUMENTS CASSETTE 10 PIECES CASSETTE	191
FRIEDMANN ROOT ELEVATOR	36	INSTRUMENTS CASSETTE 15, SURGERY KIT	191
FRIEDMANNRICHARD PERIODONTAL INSTRUMENTS	122	INSTRUMENTS CASSETTE	191
FUZZI COMPOSITE MODELLING INSTRUMENTS	85	INSTRUMENTS CASSETTE 20 DETACHABLE	191
<b>G</b>		INSTRUMENTS CASSETTE 20 DETACHABLE LIDS	191
GERALD TWEEZERS	52	INSTRUMENTS CASSETTE 20 HINGED	191
GILLES TWEEZERS	52	INSTRUMENTS CASSETTE 5 DETACHABLE LIDS	191
GILLIES NEEDLE HOLDER WITH SCISSORS	107	INSTRUMENTS CASSETTE 5 HINGED LIDS WITH BUTTON LOCK	191
GLOVES BOX	184	INSTRUMENTS CASSETTE 5 PIECES CASSETTE	191
GOLD COLOR SET DENTAL ORAL IMPLANTOLOGY SINUS	156	INSTRUMENTS CASSETTE 7 HINGED	191
SET OF 10 PCS WITH CASSETTE	174	INSTRUMENTS CASSETTE 7 PIECES CASSETTE	191
GOLD MILL (LARGE)	174	INSTRUMENTS TRAY PERFORATED	183, 191
GOLD MILL (SMALL)	174	INSTRUMENTS TRAY PERFORATED BASE	191
GOLDMAN FOX SUPER CUT SCISSORS	108	INSTRUMENTS TRAY SOLID BASE	191
GOLDMAN FOX EXPLORERS	55	INSTRUMENTS TRAYS	183
GOLDMAN FOX SCISSORS	108	INSTRUMENTS TRAYS WITH COVER	184
GOLDMAN-FOX GUM SCISSORS	104, 105	INSTRUMENTS TRAYS WITH LID	183
GOLDMAN-FOX PERIODONTAL INSTRUMENTS	121	INSTRUMENTS TRAYS WITHOUT COVER	184
GOLDMANN-FOX SCALER	64	INTRALIGAMENTAL ANESTHESIA SYRINGE	27
GRACEY AFTER FIVE	60	IRIS GUM SCISSORS	104
GRACEY CURETTES SCALER	61	IRIS SUPER CUT SCISSORS	108
GRACEY CURETTES	59, 63	IRIS GUM SCISSORS	104, 105, 108
GRACEY MINI FIVE	60	IVORY SIGVELAND RETAINERS	115
GRACEY SCALER	61	IVORY LIGHT WEIGHT RUBBER DAM CLAMPS FORCEPS	111
GRAFT MATERIAL SPOON CONCAVE TIP FOR LARGER GRAFTS	165	IVORY MATRIX RETAINERS	115
GRAFT PACKER	166	IVORY RETRACTORS	101
GRAFT PACKER, CONCAVE	166	IVORY RUBBER DAM CLAMPS FORCEPS	111
GRIPPER CARBIDE TIPPED TISSUE FORCEPS TC	169	IVANSON SPRING CALIPER FOR METAL AND PORCELAIN	126
GRITMAN CEMENT SPATULAS	95	IVANSON SPRING CALIPER FOR WAX	126
GROSS MAIER DRESSING FORCEPS	110	<b>J</b>	
GULDENER SCALER	64	JAUQUETTE SCALER	63
GUM SCISSORS	104, 108	JARABAK PILERS	143
<b>H</b>		JECT ULTIMA DENTAL SYRINGES	27
HAEMOSTATIC FORCEPS & TOWEL CLAMPS	110	JENNINGS RETRACTORS - MOUTH GAGS	100
HALSEY NEEDLE HOLDERS	106, 107	JOHNSON CONTOURING PLIERS	133
HALSTEAD MOSQUITO HAEMOSTATIC FORCEPS	110	JOHNSON PLIERS FOR CROWN BAND CONTOURING	138
HALSTEAD MOSQUITO NEEDLE HOLDERS	107	JONES TOWEL CLAMPS	110
HALTEPINZETTE WITH LOCK TWEEZERS	170	JOSEPH GUM SCISSORS	105
<b>K</b>		<b>K</b>	
		KALT NEEDLE HOLDERS TC	109
		KELLS EXTRACTING FORCEPS (AMERICAN PATTERN)	15



DESCRIPTION	PAGE#
KELLY GUM SCISSORS	104,105
KELLY HAEMOSTATIC FORCEPS	110
KELLY-RANKIN HAEMOSTATIC FORCEPS	110
KEY RINGS	94
KEYS FOR IMPLANT SCREW	177
KILNER (RAGNELL) GUM SCISSORS	105
KIRKLAND SCALER	64
KLEIN EXTRACTING FORCEPS	08
KLEINERT-KUTZ BONE ROUNDING	25
KOPP ROOTS ELEVATOR	38
KRAMER NEVINS SINUS LIFT	162
KRAMER-NEVINS PERIODONTAL INSTRUMENTS	121
KRAMER-NEVINS BONE CHISEL	119
KUEHNE TWEEZERS	170
<b>L</b>	
LA GRANGE NEEDLE HOLDERS	109
LA GRANGE GUM SCISSORS	104,105
LADMORE FILLING INSTRUMENTS	88
LADMORE IMPLANTOLOGY	149
LAGRANGE SCISSORS	108
LAGRANGE SUPER CUT SCISSORS	108
LANGENBECK RETRACTORS	103
LANGER SCALER	63
LAP & L-KEY JOINT PLIERS	139
LARGE FLUID-CUT DISTAL AND CUTTING PLIERS	133
LAXUTOR ROOT ELEVATORS	30,31
LAXUTOR ROOT ELEVATORS SET OF 9 PCS	38
LE CRON CEMENT SPATULAS	95
LECLUSE BAYONET ROOT ELEVATOR	41
LESSMANN WAX KNIVES	96
LICHTENBERG-RYDER T.C NEEDLE HOLDERS	109
LIGATURE FORMING PLIERS	134
LIGATURE INSTRUMENT SILICON GRIP	129
LIGHT WORE BIRD BEAK PLIERS WITH GROOVE	143
LINDO LEVIAN ROOT ELEVATORS	35
LINGUAL / PALATAL ARCH FORMING PLIERS	134
LINGUAL ARCH FORMING PLIERS	143
LINGUAL DENTAL PHOTOGRAPHIC MIRRORS	98
LINGUAL/PALATAL SURFACES SCALER	61
LINGUAL-LARGE DENTAL PHOTOGRAPHIC MIRRORS	98
LIP & CHEEK RETRACTOR	100
L-KEY JOINT PLIERS	139to138
LOCKLIN GUM SCISSORS	105
LONDON COLLEGE TWEEZERS	52
LONDON HOSPITAL	34
LONG IMPLANT DRILLS	176
LOOP FORMING IN FOUR LENGTH	138
LORNAHAEMOSTATIC FORCEPS	110
LUCAS IMPLANTOLOGY	168
LUER BONE ROUNDING	25
LUKS	58
LUXATING ELEVATOR / PERIOTOME	36
LUXATING ELEVATOR WITH PERIOTOME TIP	36
LUXATOR ROOT ELEVATORS	29
LUXATOR/ELEVATOR 8 PCS	191
LUXATORS SET OF 7 PCS	46
LUXATORS SET OF 7 PCS	46
LUXATORS SET OF 8 PCS	46
LUXIER ROOTS ELEVATOR	40
LUXIVATORS CURVED ATRAUMATIC LUXATION OF THE TOOTH	48
LUXIVATORS STRAIGHT ATRAUMATIC LUXATION OF THE TOOTH	48
<b>M</b>	
MAGNIFYING MOUTH MIRRORS	49
MALLET IMPLANTOLOGY	174
MALLETS - PLASTER SHEARS - SAWS & BLADES	124
MANUAL ASPIRATING SYRINGE	27
MANUAL CROWN REMOVERS	114
MARGIN TRIMMERS	66
MARGINAL RIDGE & EMBRASURE SHAPING INSTRUMENT	148
MARKLEY FILLING INSTRUMENTS	93
MARQUETTE CONDENSER	87
MARQUETTE PLUGGERS	87
MARYAM TWEEZERS	53
MATHIEU NEEDLE HOLDERS	107
MATRIX RETAINERS	115
MAXI HOLLOWARE INSTRUMENTS	187
MAYO GUM SCISSORS	104
MAYO NEEDLE HOLDERS	106
MAYO TRAY PERFORATED	183
MAYO-HEGAR NEEDLE HOLDERS	106,107
MCCALL SCALER	63
MEAD BONE ROUNDING	25
MEAD NON AUTO CLAVABLE PLASTIC TIPS	119
MEASURING LOOP	128
MEASURING LOOP CASTROVIEJO	128
MEASURING PROBE	145
MEISSNER EXTRACTING FORCEPS	26
MEMBRANE IMPLANTOLOGY	161
MEMBRANE 3POSITIONING IMPLANTOLOGY	150
MEMBRANE DETACH ELEVATOR	164
MEMBRANE ELEVATOR	164
MEMBRANE TWEEZERS	170
MERRIAM TWEEZERS	52
MERRITT & EXPLORERS	55
MERSON BAND PUSHER	129
MESH CASSETTES 5 PIECES, 7 PIECES, 10 PIECES, 15PIECES & 20 PIECES	190
MESH TRAY	190
MESH TRAYS 05 PCS	190
MESH TRAYS 07 PCS	190
MESH TRAYS 10 PCS	190
MESH TRAYS	190
MESH TRAYS 20 PCS	190
MESIAL SCALER	61
METZENBAUM GUM SCISSORS	104,106
MICOR ADSON 1X2 TISSUE FORCEPS	189
MICRO BOLEY GAUGE STRAIGHT	146
MICRO BOLEY GAUGE TWISTED	146
MICRO IRIS GUM SCISSORS	104
MICRO SURGICAL SCALPEL HANDLE TITANIUM	123
MICRO-FRIEDMAN BONE ROUNDING	25
MICRO-MINI PIN & LIGATURE CUTTER	142
MIDDELBOERF RETRACTORS	103
MILLER BONE FILES	119
MILLER ROOT ELEVATORS	37

DESCRIPTION	PAGE#
MILLER IMPLANTOLOGY	168
MILLER ROOT PLANNING FILES	119
MILLER TWEEZERS	53
MINI DISTAL AMALGAM CARRIERS	117
MINI GAUGE	127
MINI LIGATURE CUTTERS	133
MINI PIN & LIGATURE CUTTER	142
MINI PIN & LIGATURE CUTTER TC	139
MINI-LANGENBECK RETRACTORS	102
MINI-SPATEL BURNISHERS	71
MINNESOTA RETRACTORS	101
MIRROR EXPLORERS	51
MIRROR HANDLE STEM WITH RULE HANDLES FOR MOUTH MIRRORS	50
MISCELLANEOUS PRODUCTS	94
MIXING SPATULA WITH SHOVEL	165
MOLT 9 PERIODONTAL INSTRUMENTS	121
MOLT AUTO CLAVABLE PLASTIC TIPS	119
MONARI COMPOSITE MODELLING INSTRUMENTS	85
MORTENSON FILLING INSTRUMENTS	80,87
MOUTH MIRROR EXPLORERS	51
MOUTH MIRRORS	49
MOUTH PROPS	94
MOVEABLE HEIGHT GAUGE	125
MSY CEMENT SPATULAS	95
MUEHLEMAN SCALER	64,65
MULTIFUNCTIONAL METER	145
MULTI-PURPOSE PLIERS	137
<b>N</b>	
NAKAWALNI DENTAL LABORATORY INSTRUMENT	97
NANCE LOOP FORMING PLIERS	132,137
NAPKIN HOLDER WITH METAL CHAIN	94
NAPKIN HOLDER WITH PLASTIC CHAIN	94
NARROW SIOVELAND RETAINERS	115
NEEDLE HOLDERS	106to109
NEEDLE STERILIZATION TUB	186
NEUMANN GUM SCISSORS	105
NEVINS EXTRACTING FORCEPS (AMERICAN PATTERN)	14,16
NORMAL HANDLE EXPLORERS	54
NOYES GUM SCISSORS	105
<b>O</b>	
OBUSE IMPLANTOLOGY	161
OCCUSAL ANATOMY INSTRUMENT	148
OH- PERIOSTEAL ELEVATORS	120
OLSEN-HEGAR NEEDLE HOLDER WITH SCISSORS	107
OMEGA LOOP FORMING PLIERS	141
ORBAN SCALER	64
ORTHO BAND PUSHER SCALER	131
ORTHO INSTRUMENT	130
ORTHO LIGATURE DIRECTOR SINGLE ENDED	131
ORTHO LIGHT BAND PUSHER SCALER	129
ORTHO LIGHT DOUBLE ENDED LIGATURE	131
ORTHO LIGHT DOUBLE ENDED LIGATURE DIRECTOR	131
ORTHO LIGHT DOUBLE ENDED LIGATURE DIRECTOR SCALER	129,131
ORTHO LIGHT DOUBLE ENDED LIGATURE SCALER	131
ORTHO LIGHT LIGATURE DIRECTOR DOUBLE ENDED	131
ORTHO LIGHT SCALER	130
ORTHO LIGHT SCALER DOUBLE ENDED	130
ORTHO LIGHT SCALER HEAVY	130
ORTHO LIGHT SCALER SINGLE ENDED	131
ORTHO RULER / CONE SOCKET	146
ORTHODONTIC INSTRUMENTS	129to131
ORTHODONTICS PLIERS STANDS	144
ORTHOLITE LIGATURE DIRECTOR SCALER	131
OSTEOMETER MEASURING FORCEPS FOR JAW	146
OSTEOTOME CONVEXES / CONCAVE / CONIQUES	171,172
OSTEOTOME CONVEXES / CONCAVE / CONIQUES SET OF 5 PIECES	172
<b>P</b>	
P.K. THOMS WAXING INSTRUMENTS KIT	89
P.K. THOMS WAXING INSTRUMENTS SOLID HANDLE	89
PALATAL-ADULT DENTAL PHOTOGRAPHIC MIRRORS	98
PALATAL-ADULT LARGE DENTAL PHOTOGRAPHIC MIRRORS	98
PALATAL-CHILD DENTAL PHOTOGRAPHIC MIRRORS	98
PALMER RUBBER DAM CLAMPS FORCEPS	111
PARMLY EXTRACTING FORCEPS (AMERICAN PATTERN)	16
PARTSCH PERIOSTEAL ELEVATORS	120
PDL ELEVATORS	47,48
PDL ELEVATORS SET OF 7 PCS	47
PDL ELEVATORS SET OF 7 PCS	47
PDL ELEVATORS SET OF 8 PCS	48
PEDIATRIC TOOTH EXTRACTING FORCEPS SET	08
PEET SPLINTER FORCEPS	109
PERFORATED WITH RIM LOCK SET OF 06 PCS	179
PERFORATED WITH RIM LOCK SET OF 10 PCS	179
PERFORATED WITH RIM LOCK SET OF 8 PCS	179
PERIODONTAL CHISELS	118
PERIODONTAL INSTRUMENTS	120to122
PERIODONTAL KNIVES	64
PERIODONTAL PROBES	58
PERIODONTAL PROBESSPREADERS	57
PERIODONTAL SURGERY	162
PERIOSTEAL ELEVATORS	120,154
PERIOTOME ROOTS ELEVATOR	43
PERRY TWEEZERS	52,53
PHYSIC EXTRACTING FORCEPS (AMERICAN PATTERN)	18
PIN & LIGATURE CUTTER	142
PIN & LIGATURE CUTTERS TC	139
PIN HOLDING PLIER TC	140
PIN-TWEEZERS	170
PLAIN MOUTH MIRRORS	49
PLANATOR	69
PLANATOR	69
PLASTER KNIFE	86
PLASTER SHEARS	124
PLASTIC FILLING INSTRUMENTS	93
PLIER STAND PLASTIC	144
PLIER STAND STAINLESS STEEL	144
PLIERS FOR ORTHODONTICS & PROSTHETICS	132to143
PLIERS STAND STAINLESS ROD (ALUMINUM ANODIZED)	144
PLUGGERS	56
PLUGGERS/CONDENSERS	87
POSITIONING 3 IMPLANTOLOGY	161



**DESCRIPTION PAGE#**

POSITIVE POSTIONER	125
POSTERIOR BAND REMOVING PLIERS - LONG	141
POSTERIOR BAND REMOVING PLIERS - SHORT	141
POTTS "T-HANDLE" ELEVATOR-LEFT	41
POTTS "T-HANDLE" ELEVATOR-RIGHT	41
PRICHARD PERIODONTAL INSTRUMENTS	121
PROBE EXPLORERS	51
PRONG 3 CLASP ADJUSTABLE PLIER TC.	140

**Q**

QUATRO DENTAL SYRINGES	27
------------------------	----

**R**

RASPATOR, PRICHARD	165
RASPATOR, RASPATORY	120
RASPATOR, WILLIGER	120
RECTANGULAR ARCH BENDING PLIERS	143
RETRACTOR SELDIN	103
RETRACTORS - MOUTH GAGS	99,100
RETRACTORS	101to103
RHOODIUM MOUTH MIRRORS	49
RIBBON RETRACTORS	101
RIDGE MAPPING DOUBLE CALIPER	146
RING 2 DENTAL SYRINGES	27
ROCHESTER-KOCHER HAEMOSTATIC FORCEPS	110
ROCHESTER-PEAN HAEMOSTATIC FORCEPS	110
ROOT CANAL PLUGGER	56
ROOT ELEVATOR	290to48
ROOT ELEVATORS SET OF 10 PCS	46
ROOT FRAGMENT FORCEPS	09,11
ROOT PICK ELEVATOR WITH GOLD FLEXIBLE POINTS	42
ROOT PLANNING FILES	119
ROOT PLANNING FILES & BONE FILES	119
ROOTS SPLINTER FORCEPS	10
ROSEKÖNIG RETRACTORS - MOUTH GAGS	100
ROUND AND CONCAVE PLIER	139
ROUND ENDOODONTIC BOXES	186
ROUND EXCAVATOR	86
ROUND TORQUE WRENCH UNIVERSAL	177
ROUTURIER EXTRACTING FORCEPS (ENGLISH PATTERN)	03
ROUTURIER EXTRACTING FORCEPS FOR LOWER MOLARS	03
RUBBER DAM CLAMPS	113
RUBBER DAM CLAMPS FORCEPS	111
RUBBER DAM CLAMPS KITS	112
RUBBER DAM FRAME	111
RUBBER DAM TEMPLATE	112
RULER	127
RULER WEDGE	127
RUSSIAN LONG TISSUE FORCEPS	170
RUSSIAN TISSUE FORCEPS	170
RYDER NEEDLE HOLDERS	106,107

**S**

SAW BLADES	124
SCALER	59to67
SCALERS TRAYS	183
SCALPEL HANDLE STRAIGHT	123
SCALPEL HANDLES	123
SCHLEMMER ROOTS ELEVATOR	37
SCHMIDLIN PERIODONTAL INSTRUMENTS	122
SCHURE BAND SEATER	130
SCHWARTZ PLIERS	139
SEPARATING ELASTIC PLIER	111
SELDIN PERIODONTAL INSTRUMENTS	122
SELDIN ROOT ELEVATORS	29,32
SELF LOCKING GROOVED	169
SELF-LOCKING TWEEZER	53
SEMKEN TAYLOR TWEEZERS	52
SENN MUELLER RETRACTOR	102
SERRATED PLUGGERS	84
SET OF 5 PIECES	171
SHEPHERD'S HOOK EXPLORER	54
SHORT IMPLANT DRILLS	176
SICKLE SCALERS	63
SIMPLEX RETRACTOR	100
SINUS ELEVATOR & MEMBRANE LIFT TIPS	164
SINUS IMPLANTOLOGY	156,161
SINUS LEFT ELEVATOR & MEMBRANE LIFT TIPS 3 SET OF 7 PCS	164
SINUS LEFT ELEVATOR SET OF 10 PCS	158
SINUS LIFT ELEVATOR	151,157
SINUS LIFT ELEVATORS SET OF 7 PCS WITH CASSETTE	157
SINUS LIFT INSTRUMENT	150to152,161
SINUSLIFT IMPLANTOLOGY	150,161
SIOVELAND RETAINERS	115
SLIM WEINGART PLIERS	141
SLIPPER TWEEZERS	128
SMAHA ELASTOMERIC LIGATURE PLACING FORCEPS	140
SMALL MEDIUM BALL BURNISHER	149
SOFT BRUSHING KIT SET OF 6 PCS	163
SOFT WIRE AND LIGATURE CUTTER	135
SOFT WIRE AND LIGATURE CUTTER, STANDARD	135
SPATULAS	73
SPATULAS ROUNDED AND SHARP	72
SPECIAL INOX	124
SPECIAL INOX WAX KNIVES	96
SPENCER GUM SCISSORS	105
SPIRIT LAMP	186
SPLINTER FORCEPS	167
SPOON AND PUSHER (TRINOVO)	10
STANDARD LIGATURE CUTTER TC.	139
STAR TORQUE WRENCH UNIVERSAL	177
STERILIZING BOXES & INSTRUMENTS TRAYS	184
STERILIZING MESH TRAY	190
STERNBERG RETRACTORS	99
STIEGLITZ ROOT FRAGMENTS	10
STIEGLITZ SPLINTER FORCEPS	10
STOCKES RUBBER DAM CLAMPS FORCEPS	111
STOP PLIERS V-BEED	137
STRAIGHT DIRECT BOND BRACKET REMOVER TC.	140
STRIP SAW	125
SUCTION CANNULA KIT WITH 3 NEEDLES	178
SUCTION NEEDLE	178
SUCTION PIPE	178
SUCTION TUBES	178
SURGICAL TWEEZERS	52

**DESCRIPTION PAGE#**

SURGICAL BALL HOOK CRIMPING PLIERS, ANGLED	137
SURGICAL BALL HOOK CRIMPING PLIERS, STRAIGHT	137
SURGICAL CHISELS	162
SUTURE PLIERS SUTURE FORCEPS	169
SWEDISH TYPE DENTAL SYRINGES	27
SYNDESMOTOMES ROOTS ELEVATOR	43

**T**

TAMPER	166
TANNER FILLING INSTRUMENTS	88
TARTAR REMOVER INSTRUMENTS SOLID HANDLE	67
TAYLOR SCALER	63
THICK EXPLORERS	54
THREE JAW PLIERS	143
TISSUE GRAFT FORCEP SERRATED TEETH	169
TISSUE MEMBRANE FORCEPS	169,170
TISSUE PUNCH	173
TISSUE PUNCH FORCEPS TO REMOVE TISSUE AFTER USING PUNCH	170
TISSUE PUNCH SET OF FIVE PIECES	173
TISSUE PUNCH STRAIGHT & CURVED	173
TISSUE PUNCH WITH PIN SET OF FOUR PIECES	173
TITANIUM LOCKING PLIER	169
TITANIUM PERIOSTEAL ELEVATOR FOR IMPLAN	122
TOFFLEMIRE JUNIOR SIOVELAND RETAINERS	115
TOFFLEMIRE SENIOR SIOVELAND RETAINERS	115
TOMAS EXTRACTING FORCEPS (AMERICAN PATTERN)	17,22
TOWNS NEEDLE HOLDERS	109
TORQUE WRENCH & RATCHET	177
TORQUING PLIERS SET	137
TOWNER-JAQUETTE SCALER	63
TREPHINE	175
TREPHINE DRILL SET OF EIGHT PIECES	175
TRINOVO HANDLE	167
TRIO DENTAL SYRINGES	27
TRIPLE-ANGLED BACK ACTION EXCAVATOR	86
TWEED / ANGEL FORMING AND TORQUING PLIERS	134
TWEED / ANGEL RECTANGULAR ARCH FORMING AND TORQUING PLIERS	137
TWEED / ANGLE RECTANGULAR ARCH FORMING & TORQUING PLIERS	133
TWEED LOOP FORMING PLIERS	132,133,143
TWEEZERS	52,53
TWEEZERS EXPLORERS	51
TWEL-ON ELASTIC MODULE INSTRUMENT	130
TYPE C DENTAL SYRINGES	27

**U**

UNIQUE RETRACTORS	102
UNIVERSAL CROWN REMOVER	114
UNIVERSAL CURETTES	63
UNIVERSAL CUT & HOLD DISTAL END CUTTER	142
UNIVERSAL DOLPHIN STYLE WIRE BENDING & CUTTING PLIERS	132
UNIVERSAL EXTRACTING FORCEPS (AMERICAN PATTERN)	22
UNIVERSAL MINI-ELEVATOR FOR FRAGMENTS	34
UNIVERSAL NEEDLE HOLDERS	106
UNIVERSITY OF MICHIGAN EXPLORERS	55
UNIVERSITY OF NORTH CAROLINA EXPLORER	55
UPPER ROOT ELEVATOR WITH HELICOIDAL TIP	31
US ARMY RETRACTORS	102
UTILITY PLIER WEINGART TC.	140

**V**

V BEND STOP	139
V BEND STOP PLIER TC.	139
VERNIER GAUGE	127
VISTA TUNNELING INSTRUMENTS SET	155

**W**

W. BARRY ROOTS ELEVATOR	41
WAHL BONE FILES	119
WAKAI HEAT CARRIERS	57
WALDSAVHS WIRE BENDING AND CUTTING PLIERS	134
WANGEN HOLDER	99
WARD FILLING INSTRUMENTS	80
WATER SYRINGE	116
WAX KNIVES	96
WAY 6 ADJUSTABLE ANGLED SCALPEL HANDLE	123
WEIDER RETRACTOR	99
WEINGART PLIER	139
WEINGART PLIER MEDIUM	132,140
WEINGART PLIERS	141
WEINGART UTILITY PLIERS	136
WEINGART-MINI UTILITY PLIERS	136
WEITLANER RETRACTORS - MOUTH GAGS	100
WERDELSTAEDT SCALER	67
WESTCOTT BURNISHERS	68
WESTCOTT FILLING INSTRUMENTS	93
WESTON CEMENT SPATULAS	95
WHITE BURNISHERS	68
WIDE SIOVELAND RETAINERS	115
WILAND SCALER	65
WILKINS/TUFTS EXPLORER	54
WILLIAMS EXPLORERS	55
WILLIAMS FILLING INSTRUMENTS	81
WINTER EXTRACTING FORCEPS (AMERICAN PATTERN)	12
WINTER ROOTS ELEVATOR	41
WIRE BENDING & LOOP FORMING PLIERS	132
WIRE BENDING AND LOOP FORMING PLIERS ANGLED (BIRD BEAK)	136
WISDOM EXTRACTING FORCEPS	03
WISDOM ROOT ELEVATOR	36
WITZEL FORCEPS WITH ANATOMICAL HANDLE SET OF 3 PCS	09
WITZEL FORCEPS WITH GOLD SET OF 3 PCS	09
WITZEL FORCEPS WITH HOLE AND SCREW PATTERN SET OF 3 PCS	09
WITZEL ROOT FRAGMENT FORCEPS	09
WOODEN HANDLE WAX KNIVES	96

**X**

X RAY FILM HANGERS	94
--------------------	----

**Y**

YOUNG TONGUE FORCEPS	110
YOUNG WIRE BENDING & LOOP FORMING PLIERS	132
YOUNGER GOOD SCALER	63

**Z**

ZAHLE CEMENT SPATULAS	95
-----------------------	----





**Manufacturers & Exporters of Quality Dental & Surgical instruments**

**Sadaye Shehar Road. SIE Sialkot 51310 Pakistan**

**City: Sialkot-Punjab**

**Postal Code : 51310 email: [azqa.instruments@gmail.com](mailto:azqa.instruments@gmail.com)**

**[Info@azqainstruments.com](mailto:Info@azqainstruments.com)**

**Contact number +92 332 9598000**

**Whatapp/ WeChat +923087708383**

**[www.azqainstruments.com](http://www.azqainstruments.com)**

**Alibaba <https://azqa.trustpass.alibaba.com/>**